



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 03/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108186

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1844-01	GP-14-030520	1	1.15	9.61	8.99	92.7	
L1844-02	GP-15-030520	2	1.16	9.85	8.15	80.4	
L1844-03	GP-16-030520	3	1.18	9.89	8.53	84.4	
L1844-04	GP-17-030520	4	1.16	9.66	8.76	89.4	
L1844-05	GP-18-030520	5	1.14	9.53	8.32	85.6	
L1844-06	GP-19-030520	6	1.15	9.98	9.33	92.6	
L1844-07	GP-20-030520	7	1.13	9.93	9.04	89.9	
L1845-01	P001-CA42-W1-0006-01	8	1.15	9.9	8.5	84.0	
L1845-02	P001-CA43-W1-0006-01	9	1.14	9.75	8.4	84.3	
L1845-03	P001-CA44-W1-0006-01	10	1.13	9.77	8.36	83.7	
L1845-04	P001-CB41-W1-0006-01	11	1.13	9.87	8.76	87.3	
L1845-05	P001-CB41-W2-0006-01	12	1.18	9.8	8.39	83.6	
L1845-06	P001-CB44-W1-0006-01	13	1.15	9.97	9.26	92.0	
L1845-07	P001-CB44-W2-0006-01	14	1.16	9.74	8.63	87.1	
L1845-08	P001-CB45-W1-0006-01	15	1.13	9.72	9.02	91.9	
L1845-09	P001-CB45-W2-0006-01	16	1.13	9.64	9.08	93.4	
L1845-10	P001-CB45-W2-0006-01D	17	1.13	9.64	9.08	93.4	
L1845-11	P001-CB45-W2-0006-01S	18	1.13	9.64	9.08	93.4	
L1845-12	P001-CB45-W2-0006-02	19	1.15	9.7	9.18	93.9	
L1845-13	P001-CC40-W1-0006-01	20	1.14	9.61	8.55	87.5	
L1845-14	P001-CC40-W2-0006-01	21	1.13	9.96	9.07	89.9	
L1845-15	P001-CC46-W1-0006-01	22	1.18	9.57	8.02	81.5	
L1845-16	P001-CC46-W2-0006-01	23	1.16	9.63	9.06	93.3	
L1845-17	P001-CC47-W1-0006-01	24	1.16	9.52	7.78	79.2	
L1845-18	P001-CD40-W1-0006-01	25	1.14	9.97	8.18	79.7	
L1845-19	P001-CD40-W2-0006-01	26	1.14	9.95	8.97	88.9	
L1845-20	P001-CD48-W1-0006-01	27	1.13	9.87	9.11	91.3	
L1845-21	P001-CD48-W2-0006-01	28	1.13	9.74	7.58	74.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 03/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108186

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1846-01	P001-CE41-W1-0006-01	29	1.15	9.59	8.82	90.9	
L1846-02	P001-CE49-W1-0006-01	30	1.15	9.84	8.88	89.0	
L1846-03	P001-CE50-W1-0006-01	31	1.13	9.68	7.49	74.4	
L1846-04	P001-CE50-W2-0006-01	32	1.13	9.97	9.2	91.3	
L1846-05	P001-CF41-W1-0006-01	33	1.12	9.82	8.12	80.5	
L1846-06	P001-CF50-W1-0006-01	34	1.14	9.85	8.17	80.7	
L1846-07	P001-CG41-W1-0006-01	35	1.12	9.98	8.78	86.5	
L1846-08	P001-CG41-W2-0006-01	36	1.12	9.82	9.03	90.9	
L1846-10	P001-CG43-W1-0006-01	37	1.15	9.51	8.61	89.2	
L1846-11	P001-CG43-W1-0006-01D	38	1.15	9.51	8.61	89.2	
L1846-12	P001-CG43-W1-0006-01S	39	1.15	9.51	8.61	89.2	
L1846-13	P001-CG43-W1-0006-02	40	1.11	9.57	8.65	89.1	
L1846-14	P001-CG44-W1-0006-01	41	1.14	9.88	8.66	86.0	
L1846-15	P001-CG45-W1-0006-01	42	1.12	9.98	9.31	92.4	
L1846-16	P001-CG50-W1-0006-01	43	1.13	9.96	7.47	71.8	
L1846-17	P001-CH45-W1-0006-01	44	1.14	9.97	8.64	84.9	
L1846-21	P001-CH50-W1-0006-02	45	1.16	9.78	9.15	92.7	
L1846-22	P001-CI45-W1-0006-01	46	1.13	9.87	8.72	86.8	
L1846-23	P001-CI50-W1-0006-01	47	1.18	9.98	7.98	77.3	
L1847-01	P001-CJ45-W1-0006-01	48	1.12	9.63	8.87	91.1	
L1847-02	P001-CJ50-W1-0006-01	49	1.16	9.97	8.29	80.9	
L1847-03	P001-CK45-W1-0006-01	50	1.16	9.75	8.95	90.7	
L1847-04	P001-CK45-W2-0006-01	51	1.13	9.52	8.92	92.8	
L1847-05	P001-CK45-W3-0006-01	52	1.13	9.85	8.66	86.4	
L1847-06	P001-CK46-W1-0006-01	53	1.15	9.66	6.95	68.2	
L1847-07	P001-CK46-W2-0006-01	54	1.16	9.8	7.48	73.1	
L1847-08	P001-CK47-W1-0006-01	55	1.15	9.82	7.57	74.0	
L1847-09	P001-CK47-W2-0006-01	56	1.15	9.67	7.82	78.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 03/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108186

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1847-10	P001-CK47-W3-0006-01	57	1.14	9.87	7.78	76.1	
L1847-11	P001-CK48-W1-0006-01	58	1.13	9.65	7.98	80.4	
L1847-12	P001-CK48-W2-0006-01	59	1.14	9.74	8.91	90.3	
L1847-13	P001-CK48-W3-0006-01	60	1.13	9.92	9.33	93.3	
L1847-14	P001-CK49-W1-0006-01	61	1.13	9.78	7.64	75.3	
L1847-15	P001-CK50-W1-0006-01	62	1.13	9.74	8.91	90.4	
L1847-16	P001-CK50-W2-0006-01	63	1.14	9.86	7.87	77.2	
L1847-17	P001-CK50-W3-0006-01	64	1.14	9.65	7.23	71.6	
L1847-18	P001-CK50-W4-0006-01	65	1.14	9.91	8.35	82.2	
L1847-19	P001-CH50-W1-0006-01	66	1.14	9.86	9.2	92.4	
L1847-20	P001-CH50-W1-0006-01D	67	1.14	9.86	9.2	92.4	
L1847-21	P001-CH50-W1-0006-01S	68	1.14	9.86	9.2	92.4	
L1849-01	VNJ-230-218-211-216-242	69	1.14	9.94	9.04	89.8	
L1849-03	VNJ-215-72-12018	70	1.13	9.71	8.4	84.7	
L1852-01	TP-1-S1	71	1.15	9.95	9.8	98.3	
L1852-02	TP-1-S2	72	1.14	9.94	9.56	95.7	
L1852-03	TP-1-S	73	1.12	9.7	9.06	92.5	
L1852-05	TP-2-S1	74	1.1	9.89	8.61	85.4	
L1852-06	TP-2-S2	75	1.16	9.69	8.42	85.1	
L1852-07	TP-2-S	76	1.17	9.68	8.35	84.4	
L1852-09	TP-2-D1	77	1.13	9.78	8.02	79.7	
L1852-10	TP-2-D2	78	1.17	9.91	8.41	82.8	
L1852-11	TP-2-D	79	1.12	9.55	7.71	78.2	
L1852-13	TP-3-S1	80	1.13	9.51	7.56	76.7	
L1852-14	TP-3-S2	81	1.1	9.6	7.49	75.2	
L1852-15	TP-3-S	82	1.17	9.77	7.77	76.7	
L1852-17	TP-3-D1	83	1.15	9.82	8.51	84.9	
L1852-18	TP-3-D2	84	1.16	9.52	8.2	84.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 03/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108186

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1852-19	TP-3-D	85	1.12	9.55	8.32	85.4	
L1852-21	TP-4-S1	86	1.17	9.95	9.69	97.0	
L1852-22	TP-4-S2	87	1.16	9.97	9.55	95.2	
L1852-23	TP-4-S	88	1.12	9.94	9.38	93.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$