



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 03/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 03/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108247

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1761-06	SU-01-031220	1	1.12	9.51	8.58	88.9	
L1761-07	SU-01-031220	2	1.14	9.83	8.97	90.1	
L1761-08	SU-01-031220-1	3	1.17	9.9	9.12	91.1	
L1761-09	SU-01-031220-2	4	1.17	9.78	8.63	86.6	
L1762-13	TR-04-031220	5	1.15	9.64	9.33	96.3	
L1762-14	TR-04-031220	6	1.16	9.98	9.36	93.0	
L1893-13	TP-4 (3-4)	7	1.18	9.78	8.31	82.9	
L1902-01	VNJ-229-13	8	1.16	9.56	8.00	81.4	
L1907-01	TP-6 (1-2)	9	1.16	9.74	8.67	87.5	
L1907-02	TP-6-SA	10	1.11	9.54	8.63	89.2	
L1907-03	TP-6-SB	11	1.15	9.63	8.65	88.4	
L1907-05	TP-9 (3-4)	12	1.17	9.93	8.48	83.4	
L1907-06	TP-9-SA	13	1.14	9.78	8.66	87.0	
L1912-01	Z1-1	14	1.18	9.98	9.16	90.7	
L1912-02	Z1-1	15	1.18	9.98	9.16	90.7	
L1912-03	Z1-1A	16	1.14	9.51	8.96	93.4	
L1912-04	Z1-1B	17	1.19	9.81	9.2	92.9	
L1912-05	Z3-S1	18	1.18	9.69	9.06	92.6	
L1912-06	Z3-S1	19	1.18	9.69	9.06	92.6	
L1912-07	Z3-S1A	20	1.18	9.58	9.04	93.6	
L1912-08	Z3-S1B	21	1.17	9.75	9.22	93.8	
L1912-09	Z3-S2	22	1.17	9.75	9.32	95.0	
L1912-10	Z3-S2	23	1.17	9.75	9.32	95.0	
L1912-11	Z3-S2-A	24	1.13	9.9	9.08	90.6	
L1912-12	Z3-S2-B	25	1.1	9.69	8.81	89.8	
L1912-13	Z1-2	26	1.15	9.96	9.42	93.9	
L1912-14	Z1-2	27	1.15	9.96	9.42	93.9	
L1912-15	Z1-2A	28	1.11	9.81	8.85	89.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 03/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 03/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108247

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1912-16	Z1-2B	29	1.18	9.96	9.00	89.1	
L1912-17	Z1-3	30	1.17	9.8	8.76	87.9	
L1912-18	Z1-3	31	1.17	9.8	8.76	87.9	
L1912-19	Z1-3A	32	1.15	9.94	8.98	89.1	
L1912-20	Z1-3B	33	1.13	9.91	9.00	89.6	
L1912-21	Z1-4	34	1.15	9.79	8.76	88.1	
L1912-22	Z1-4	35	1.15	9.79	8.76	88.1	
L1912-23	Z1-4A	36	1.17	9.94	8.81	87.1	
L1912-24	Z1-4B	37	1.18	9.97	8.66	85.1	
L1912-25	Z1-5	38	1.11	9.94	8.35	82.0	
L1912-26	Z1-5	39	1.11	9.94	8.35	82.0	
L1912-27	Z1-5A	40	1.14	9.93	9.02	89.6	
L1912-28	Z1-5B	41	1.17	9.97	8.85	87.3	
L1912-29	Z1-6	42	1.11	9.76	9.06	91.9	
L1912-30	Z1-6	43	1.11	9.76	9.06	91.9	
L1912-31	Z1-6A	44	1.1	9.81	9.1	91.8	
L1912-32	Z1-6B	45	1.18	9.59	8.91	91.9	
L1913-01	MK-031220-0-1-GRAB-B1	46	1.12	9.66	8.39	85.1	
L1913-02	MK-031220-0-2-COMP-B1	47	1.15	9.96	8.89	87.9	
L1913-03	MK-031220-2-5-COMP-B1	48	1.12	9.5	8.44	87.4	
L1913-04	MK-031220-4-5-GRAB-B1	49	1.12	9.52	8.28	85.2	
L1913-05	MK-031220-0-1-GRAB-B3	50	1.17	9.68	8.99	91.9	
L1913-06	MK-031220-0-2-COMP-B3	51	1.17	9.54	8.63	89.1	
L1913-07	MK-031220-2-5-COMP-B3	52	1.14	9.59	9.11	94.3	
L1913-08	MK-031220-4-5-GRAB-B3	53	1.18	9.54	9.08	94.5	
L1913-09	MK-031220-0-1-GRAB-B2	54	1.12	9.89	8.65	85.9	
L1913-10	MK-031220-0-2-COMP-B2	55	1.18	9.57	8.46	86.8	
L1913-11	MK-031220-2-5-COMP-B2	56	1.14	9.93	9.18	91.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 03/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 03/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108247

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L1913-12	MK-031220-4-5-GRAB-B2	57	1.18	9.92	9.18	91.5	
L1913-13	MK-031220-0-1-GRAB-B4	58	1.16	9.64	8.84	90.6	
L1913-14	MK-031220-0-2-COMP-B4	59	1.11	9.77	9.03	91.5	
L1913-15	MK-031220-2-5-COMP-B4	60	1.17	9.89	8.84	88.0	
L1913-16	MK-031220-4-5-GRAB-B4	61	1.16	9.54	8.5	87.6	
L1913-17	MK-031220-GRAB-DP	62	1.13	9.87	9.15	91.8	
L1913-18	MK-031220-COMP-DP	63	1.13	9.98	9.11	90.2	
L1914-01	32292-32293	64	1.00	2.00	2.00	100.0	OILY DEBRIS, 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$