



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 3/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 03/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 09:00
Out Date: 03/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108283

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L1938-01	CB721	1	1.17	9.59	9.51	99.0	
L1938-02	CB721MS	2	1.17	9.59	9.51	99.0	
L1938-03	CB721MSD	3	1.17	9.59	9.51	99.0	
L1938-04	CB722	4	1.15	9.74	9.68	99.3	
L1938-05	CB723	5	1.16	9.6	9.5	98.8	
L1938-06	CB724	6	1.12	9.53	9.42	98.7	
L1938-07	CB725	7	1.17	9.92	9.8	98.6	
L1938-08	CB726	8	1.12	9.66	9.55	98.7	
L1938-09	CB727	9	1.13	9.68	9.63	99.4	
L1938-10	CB728	10	1.15	9.85	9.78	99.2	
L1938-11	CB778	11	1.13	9.82	9.7	98.6	
L1938-12	CB779	12	1.19	9.62	9.5	98.6	
L1938-13	CB780	13	1.13	9.56	9.53	99.6	
L1938-14	CB788	14	1.18	9.77	9.64	98.5	
L1938-15	CB789	15	1.19	9.88	9.75	98.5	
L1938-16	CB790	16	1.17	9.7	9.6	98.8	
L1938-17	CB791	17	1.16	9.64	9.51	98.5	
L1938-18	CB792	18	1.17	9.53	9.4	98.4	
L1938-19	CB793	19	1.1	9.76	9.64	98.6	
L1938-20	CB794	20	1.14	9.58	9.45	98.5	
L1938-21	CB795	21	1.13	9.58	9.43	98.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$