



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108520

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2074-01	S1-01-A	1	1.14	9.52	6.12	59.4	
L2075-01	S2-01-A	2	1.15	9.71	8.00	80.0	
L2075-02	S2-02-A	3	1.16	9.76	8.05	80.1	
L2075-03	S2-03-A	4	1.14	9.81	8.62	86.3	
L2075-04	S2-04-A	5	1.16	9.83	8.49	84.5	
L2075-05	S2-05-A	6	1.17	9.88	7.44	72.0	
L2075-06	S2-06-A	7	1.15	9.56	8.43	86.6	
L2075-19	DUP-1	8	1.18	9.69	7.98	79.9	
L2076-01	S3-01-A	9	1.15	9.94	8.42	82.7	
L2076-02	S3-02-A	10	1.14	9.6	8.26	84.2	
L2076-03	S3-03-A	11	1.16	9.71	8.12	81.4	
L2076-04	S3-04-A	12	1.18	9.77	7.8	77.1	
L2076-05	S3-05-A	13	1.14	9.65	8.23	83.3	
L2076-06	S3-06-A	14	1.12	9.78	7.96	79.0	
L2076-20	DUP-4	15	1.12	9.53	8.18	83.9	
L2077-01	S4-01-A	16	1.15	9.95	8.47	83.2	
L2077-02	S4-02-A	17	1.15	9.56	8.06	82.2	
L2077-03	S4-05-A	18	1.18	9.68	7.76	77.4	
L2077-04	S4-06-A	19	1.13	9.58	7.94	80.6	
L2078-01	S5-01-A	20	1.11	9.66	8.77	89.6	
L2078-02	S5-02-A	21	1.13	9.52	8.17	83.9	
L2079-01	S6-01-A	22	1.13	9.77	8.23	82.2	
L2079-02	S6-02-A	23	1.1	9.91	7.79	75.9	
L2079-03	S6-03-A	24	1.18	9.97	9.04	89.4	
L2079-04	S6-04-A	25	1.18	9.72	7.86	78.2	
L2079-05	S6-05-A	26	1.19	9.8	8.38	83.5	
L2079-16	DUP-7	27	1.13	9.59	8.01	81.3	
L2080-02	S7-02-A	28	1.17	9.58	7.69	77.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108520

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2080-03	S7-03-A	29	1.14	9.95	8.38	82.2	
L2080-10	S7-01-A	30	1.18	9.5	6.83	67.9	
L2081-01	S8-01-A	31	1.16	9.7	8.26	83.1	
L2081-02	S8-02-A	32	1.18	9.68	8.9	90.8	
L2081-03	S8-03-A	33	1.18	9.78	8.42	84.2	
L2081-04	S8-04-A	34	1.11	9.72	8.00	80.0	
L2081-05	S8-05-A	35	1.14	9.93	8.73	86.3	
L2081-06	S8-06-A	36	1.14	9.83	7.73	75.8	
L2081-07	S8-07-A	37	1.17	9.9	8.78	87.2	
L2081-08	S8-08-A	38	1.17	9.69	9.14	93.5	
L2081-09	S8-09-A	39	1.13	9.87	9.1	91.2	
L2081-10	S8-10-A	40	1.18	9.75	8.32	83.3	
L2081-26	DUP-10	41	1.14	9.84	8.77	87.7	
L2082-01	S9-01-A	42	1.14	9.96	8.3	81.2	
L2082-02	S9-02-A	43	1.15	9.94	7.67	74.2	
L2083-01	S10-01-A	44	1.17	9.92	5.61	50.7	
L2083-02	S10-02-A	45	1.13	9.63	7.55	75.5	
L2083-03	S10-03-A	46	1.14	9.9	6.19	57.6	
L2083-04	S10-04-A	47	1.11	9.82	7.8	76.8	
L2083-05	S10-05-A	48	1.17	9.98	8.19	79.7	
L2083-06	S10-06-A	49	1.12	9.97	8.22	80.2	
L2083-07	S10-07-A	50	1.14	9.56	7.58	76.5	
L2083-15	DUP-17	51	1.13	9.83	8.14	80.6	
L2084-01	S11-01-A	52	1.15	9.96	8.63	84.9	
L2084-02	S11-02-A	53	1.11	9.96	8.48	83.3	
L2084-03	S11-03-A	54	1.17	9.87	8.6	85.4	
L2085-01	S12-01-A	55	1.16	9.8	8.52	85.2	
L2085-02	S12-02-A	56	1.17	9.94	8.25	80.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108520

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2085-03	S12-03-A	57	1.1	9.71	7.68	76.4	
L2085-04	S12-04-A	58	1.13	9.66	7.92	79.6	
L2085-05	S12-05-A	59	1.17	9.75	7.28	71.2	
L2085-06	S12-06-A	60	1.15	9.51	6.36	62.3	
L2085-07	S12-07-A	61	1.13	9.73	6.9	67.1	
L2085-08	S12-08-A	62	1.13	9.61	6.33	61.3	
L2085-09	S12-09-A	63	1.16	9.96	6.87	64.9	
L2085-10	S12-10-A	64	1.13	9.94	7.73	74.9	
L2085-20	DUP-14	65	1.13	9.5	7.98	81.8	
L2086-01	S13-01-A	66	1.15	9.96	7.75	74.9	
L2086-02	S13-02-A	67	1.16	9.82	8.85	88.8	
L2086-03	S13-03-A	68	1.18	9.98	8.79	86.5	
L2087-01	S14-01-A	69	1.13	9.59	7.7	77.7	
L2087-02	S14-02-A	70	1.17	9.69	8.09	81.2	
L2087-07	DUP-15	71	1.19	9.97	8.58	84.2	
L2089-01	SB-1D- (0-2)	72	1.15	9.74	7.93	78.9	
L2089-02	SB-2D- (0-2)	73	1.18	9.88	8.04	78.9	
L2089-03	SB-3J- (0-2)	74	1.17	9.92	8.02	78.3	
L2089-04	SB-3J- (5-7)	75	1.18	9.97	8.08	78.5	
L2089-05	SB-4I- (0-2)	76	1.12	9.85	8.06	79.5	
L2096-01	RO-1-2-TP	77	1.16	9.75	8.68	87.5	
L2096-02	RO-1-2-TP	78	1.16	9.75	8.68	87.5	
L2096-03	VNJ-219	79	1.14	9.69	8.83	89.9	
L2097-17	40524	80	1.00	2.00	2.00	100.0	OILY DEBRIS, 100% SOLIDS
L2097-20	COMP-7	81	1.13	9.74	9.09	92.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$