



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108536

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2100-01	239-SB-03-5-6	1	1.15	9.96	8.47	83.1	
L2100-02	239-SB-03-10-11	2	1.17	9.75	8.76	88.5	
L2100-03	239-SB-04-5-6	3	1.17	9.98	8.67	85.1	
L2100-04	239-SB-04-10-11	4	1.17	9.94	9.27	92.4	
L2100-05	239-01-0-21	5	1.12	9.61	8.46	86.5	
L2100-06	239-05060708-0-6	6	1.17	9.95	8.7	85.8	
L2100-07	239-105060708-0-6	7	1.18	9.96	8.85	87.4	
L2100-08	239-020304-0-11	8	1.12	9.63	8.29	84.3	
L2101-01	239-WC-S-SB-01	9	1.11	9.82	8.63	86.3	
L2101-02	239-WC-S-SB-BC	10	1.15	9.94	8.69	85.8	
L2101-03	239-WC-S-SB-020304	11	1.17	9.62	8.43	85.9	
L2102-01	RT-46-MEDIAN	12	1.15	9.58	8.84	91.2	
L2102-02	RT-46-MEDIAN	13	1.15	9.58	8.84	91.2	
L2103-01	TP-60-61	14	1.14	6.3	5.44	83.3	
L2103-02	TP-61-2.0-2.5	15	1.13	9.9	8.59	85.1	
L2103-03	TP-61-1.0-1.5	16	1.12	9.58	8.21	83.8	
L2103-04	TP-61-3.5-4.0	17	1.15	9.95	8.44	82.8	
L2103-05	TP-60-61-COMP	18	1.12	9.91	8.08	79.2	
L2108-01	TP-70-71	34	1.17	6.2	5.3	82.1	
L2108-02	TP-71-2.0-2.5	35	1.18	9.65	8.06	81.2	
L2108-03	TP-71-7.0-7.5	36	1.14	9.82	8.3	82.5	
L2108-04	TP-71-12.5-13.0	37	1.13	9.7	8.08	81.1	
L2108-05	TP-70-71-COMP	38	1.19	9.79	8.22	81.7	
L2122-01	SU-01-040120	39	1.15	9.96	9.04	89.6	
L2122-02	SU-01-040120	40	1.13	9.89	8.56	84.8	
L2122-03	SU-01-040120-1	41	1.17	9.97	9.3	92.4	
L2122-04	SU-01-040120-2	42	1.13	9.94	8.28	81.2	
L2125-01	A-1	19	1.13	9.71	9.14	93.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/2/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 04/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108536

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2125-02	A-1A	20	1.17	9.71	9.18	93.8	
L2127-01	S-1	21	1.14	9.96	9.22	91.6	
L2127-02	S-2	22	1.11	9.5	8.71	90.6	
L2127-03	S-3	23	1.13	9.53	8.68	89.9	
L2127-04	S-4	24	1.17	9.75	9.00	91.3	
L2127-05	S-5	25	1.16	9.95	9.2	91.5	
L2127-06	S-6	26	1.17	9.97	9.14	90.6	
L2127-07	S-7	27	1.18	9.81	9.04	91.1	
L2127-08	S-8	28	1.18	9.92	9.07	90.3	
L2127-09	S-9	29	1.17	9.78	9.03	91.3	
L2127-10	S-10	30	1.17	9.83	8.98	90.2	
L2127-11	S-11	31	1.15	9.52	8.74	90.7	
L2130-01	400318075	32	1.18	9.98	1.52	3.9	SLUDGE SAMPLE
L2131-01	400320036	33	1.19	9.77	1.37	2.1	SLUDGE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$