



PERCENT SOLID

Supervisor: Daniel
Analyst: JIGNESH
Date: 4/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 04/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 04/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108571

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2110-01	AU-05-040220	1	1.12	9.61	9.04	93.3	
L2110-02	AU-05-040220	2	1.15	9.96	9.00	89.1	
L2110-03	AU-05-040220	3	1.14	9.63	8.83	90.6	
L2153-01	VNJ-231	28	1.15	9.52	8.66	89.7	
L2153-02	VNJ-231	29	1.18	9.88	9.17	91.8	
L2153-03	72-11991	30	1.15	9.56	9.24	96.2	
L2153-04	72-11991	31	1.13	9.77	9.2	93.4	
L2153-05	RBR-2520699	32	1.13	9.56	8.39	86.1	
L2153-06	RBR-2520699	33	1.18	9.77	8.62	86.6	
L2154-01	VNJ-212	34	1.18	9.59	8.16	83.0	
L2157-01	VNJ-226	35	1.14	9.56	3.07	22.9	
L2157-02	VNJ-226	36	1.13	9.9	3.16	23.1	
L2162-01	TP-54-53	4	1.14	5.98	5.51	90.3	
L2162-03	TP-53-2.5-3.0	5	1.18	9.72	8.97	91.2	
L2162-04	TP-53-1.5-2.0	6	1.13	9.77	8.97	90.7	
L2162-05	TP-54-53-COMP	7	1.14	9.94	9.16	91.1	
L2162-08	TP-50	8	1.15	7.38	6.22	81.4	
L2162-09	TP-50-4.5-5.0	9	1.16	9.97	8.2	79.9	
L2162-10	TP-50-3.5-4.0	10	1.17	9.92	8.77	86.9	
L2162-11	TP-50-2.5-3.0	11	1.18	9.95	8.81	87.0	
L2162-12	TP-50-COMP	12	1.18	9.92	8.64	85.4	
L2164-01	ROLL-OFF-GRAB	37	1.13	9.59	8.46	86.6	
L2164-02	ROLL-OFF-COMP	38	1.15	9.94	8.9	88.2	
L2165-01	2ND-FLOOR	39	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100% SOLIDS
L2165-03	3RD-FLOOR	40	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100% SOLIDS
L2166-01	TP-47	13	1.13	7.24	6.35	85.4	
L2166-02	TP-47-2.5-3.0	14	1.15	9.94	9.11	90.6	
L2166-03	TP-47-2.0-2.5	15	1.14	9.93	8.73	86.3	



PERCENT SOLID

Supervisor: Daniel
Analyst: JIGNESH
Date: 4/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 04/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 04/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108571

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2166-04	TP-47-1.5-2.0	16	1.15	9.5	8.09	83.1	
L2166-05	TP-47-COMP	17	1.14	9.55	8.36	85.9	
L2166-08	TP-51-1.0-1.5	18	1.14	9.96	9.02	89.3	
L2166-09	TP-49-2.0-2.5	19	1.19	9.79	8.93	90.0	
L2166-10	TP-48-2.0-2.5	20	1.15	9.96	9.41	93.8	
L2166-11	TP-46	21	1.17	5.82	4.84	78.9	
L2166-12	TP-52-46-COMP	22	1.15	9.94	9.31	92.8	
L2166-15	TP-45	23	1.14	7.37	7.2	97.3	
L2166-16	TP-45-3.5-4.0	24	1.12	9.92	9.21	91.9	
L2166-17	TP-45-3.0-3.5	25	1.18	9.67	8.36	84.6	
L2166-18	TP-45-2.5-3.0	26	1.15	9.96	8.7	85.7	
L2166-19	TP-45-COMP	27	1.18	9.71	8.54	86.3	
L2167-15	SVOC-GPC-BLANK	41	1.00	2.00	2.00	100.0	
L2167-16	PEST-GPC-BLANK	42	1.00	2.00	2.00	100.0	
L2167-17	PEST-GPC-BLANK-SPIKE	43	1.00	2.00	2.00	100.0	
L2167-18	PCB-GPC-BLANK	44	1.00	2.00	2.00	100.0	
L2167-19	PCB-GPC-BLANK-SPIKE	45	1.00	2.00	2.00	100.0	
L2167-20	SVOC-GCP2-BLANK	46	1.00	2.00	2.00	100.0	
L2167-21	PEST-GCP2-BLANK	47	1.00	2.00	2.00	100.0	
L2167-22	PEST-GCP2-BLANK-SPIKE	48	1.00	2.00	2.00	100.0	
L2167-23	PCB-GCP2-BLANK	49	1.00	2.00	2.00	100.0	
L2167-24	PCB-GCP2-BLANK-SPIKE	50	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$