



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:50
In Date: 04/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:32
Out Date: 04/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108597

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2156-01	MDBF99	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
L2156-02	MDBFA0	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
L2156-03	MDBFA1	3	1.15	9.76	8.29	82.9	
L2156-04	MDBFA2	4	1.15	9.96	8.54	83.9	
L2156-05	MDBFA3	5	1.18	9.69	7.95	79.6	
L2156-06	MDBFA3D	6	1.18	9.69	7.95	79.6	
L2156-07	MDBFA3S	7	1.18	9.69	7.95	79.6	
L2156-08	MDBFA4	8	1.15	9.96	8.66	85.2	
L2156-09	MDBFA5	9	1.14	9.86	8.18	80.7	
L2156-10	MDBFA6	10	1.12	9.92	8.01	78.3	
L2156-11	MDBFA7	11	1.16	9.97	8.26	80.6	
L2156-12	MDBFA8	12	1.18	9.58	8.19	83.5	
L2156-13	MDBFA9	13	1.18	9.68	8.32	84.0	
L2156-14	MDBFB0	14	1.14	9.65	8.27	83.8	
L2156-15	MDBFB1	15	1.13	9.87	8.18	80.7	
L2156-16	MDBFB2	16	1.14	9.9	8.39	82.8	
L2156-17	MDBFB3	17	1.19	9.93	7.19	68.6	
L2156-18	MDBFB4	18	1.19	9.93	8.17	79.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$