



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 4/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:04
Out Date: 04/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108631

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2180-05	NEWPORT-TP-7FT-B	1	1.13	9.66	8.3	84.1	
L2185-02	30941	2	1.13	9.76	9.72	99.5	
L2185-03	10:1-COMP	3	1.13	9.51	8.01	82.1	
L2186-01	COMP-8:1-SOLIDS	4	1.00	2.00	2.00	100.0	OILY DEBRIS, 100% SOLIDS
L2196-01	FE-SB-32-(8.5-9)	5	1.16	9.68	8.91	91.0	
L2196-02	FE-SB-35-(7-7.5)	6	1.14	9.9	8.02	78.5	
L2196-03	FE-SB-29-(5.5-6)	7	1.14	9.66	8.22	83.1	
L2196-04	SB-30-(8-8.5)	8	1.12	9.52	7.71	78.5	
L2196-05	SB-31-(5.5-6)	9	1.13	9.8	7.68	75.5	
L2199-01	TP-39-38	10	1.12	6.36	5.58	85.1	
L2199-02	TP-39-38-4.5-5.0	11	1.12	9.62	8.53	87.2	
L2199-03	TP-39-38-3.5-4.0	12	1.15	9.96	8.62	84.8	
L2199-04	TP-39-38-2.5-3.0	13	1.13	9.53	8.2	84.2	
L2199-05	TP-39-38-COMP	14	1.11	9.7	8.37	84.5	
L2199-08	TP-37-36	15	1.12	5.56	5.05	88.5	
L2199-09	TP-37-36-2.5-3.0	16	1.18	9.56	8.57	88.2	
L2199-10	TP-37-36-2.0-2.5	17	1.18	9.71	8.55	86.4	
L2199-11	TP-37-36-1.5-2.0	18	1.18	9.55	8.04	82.0	
L2199-12	TP-37-36-COMP	19	1.13	9.58	8.48	87.0	
L2201-01	CG-SUB	20	1.16	9.86	8.47	84.0	
L2201-02	CG-SUB-A	21	1.15	9.97	8.58	84.2	
L2201-03	CG-SUB-B	22	1.17	9.77	8.43	84.4	
L2201-04	CG-SUB-VO	23	1.14	9.51	8.34	86.0	
L2202-01	400402720.1	24	1.11	9.68	2.92	21.1	SLUDGE SAMPLE
L2203-01	400401739.0	25	1.15	9.96	4.87	42.2	SLUDGE SAMPLE
L2204-01	400406711.0	26	1.16	9.98	3.07	21.7	SLUDGE SAMPLE
L2205-01	JC-02-04082020-A	27	1.15	9.94	9.19	91.5	
L2205-02	JC-02-04082020-A	28	1.13	9.88	9.15	91.7	



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 4/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:04
Out Date: 04/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108631

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2205-03	JC-02-04082020-1	29	1.13	9.68	9.00	92.0	
L2205-04	JC-02-04082020-2	30	1.11	9.74	8.9	90.3	
L2205-05	JC-02-04082020-3	31	1.12	9.64	9.23	95.2	
L2205-06	JC-02-04082020-4	32	1.12	9.7	8.63	87.5	
L2205-07	JC-02-04082020-5	33	1.13	9.64	8.66	88.5	
L2205-09	JC-02-04082020-B	34	1.13	9.68	8.89	90.8	
L2205-10	JC-02-04082020-B	35	1.13	9.5	8.92	93.1	
L2205-11	JC-02-04082020-6	36	1.12	9.97	9.16	90.8	
L2205-12	JC-02-04082020-7	37	1.14	9.84	9.03	90.7	
L2205-13	JC-02-04082020-8	38	1.19	9.75	9.24	94.0	
L2205-14	JC-02-04082020-9	39	1.1	9.88	8.75	87.1	
L2205-15	JC-02-04082020-10	40	1.13	9.92	8.88	88.2	
L2205-17	JC-02-04082020-C	41	1.13	9.67	8.82	90.0	
L2205-18	JC-02-04082020-C	42	1.16	9.88	9.14	91.5	
L2205-19	JC-02-04082020-11	43	1.18	9.89	9.13	91.3	
L2205-20	JC-02-04082020-12	44	1.13	9.6	8.9	91.7	
L2205-21	JC-02-04082020-13	45	1.11	9.72	8.72	88.4	
L2205-22	JC-02-04082020-14	46	1.13	9.79	8.67	87.1	
L2208-01	T-2	47	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$