



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 4/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:50
In Date: 04/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108633

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2191-01	CB796	1	1.15	9.75	9.65	98.8	
L2191-02	CB796MS	2	1.15	9.75	9.65	98.8	
L2191-03	CB796MSD	3	1.15	9.75	9.65	98.8	
L2191-04	CB797	4	1.12	9.58	9.48	98.8	
L2191-05	CB798	5	1.15	9.68	9.58	98.8	
L2191-06	CB799	6	1.18	9.66	9.57	98.9	
L2191-07	CB7A0	7	1.16	9.6	9.54	99.3	
L2191-08	CB7A1	8	1.13	9.58	9.52	99.3	
L2191-09	CB7A2	9	1.16	9.75	9.6	98.3	
L2191-10	CB7A3	10	1.18	9.72	9.57	98.2	
L2191-11	CB7A4	11	1.12	9.58	9.42	98.1	
L2191-12	CB7A5	12	1.18	9.57	9.43	98.3	
L2191-13	CB7A6	13	1.19	9.65	9.56	98.9	
L2191-14	CB7A7	14	1.15	9.9	9.81	99.0	
L2191-15	CB7A8	15	1.13	9.81	9.8	99.9	
L2191-16	CB7B6	16	1.17	9.74	9.6	98.4	
L2191-17	CB7B7	17	1.14	9.5	9.38	98.6	
L2191-18	CB7B8	18	1.18	9.58	9.48	98.8	
L2191-19	CB7B9	19	1.17	9.64	9.5	98.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$