



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 4/9/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 04/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108634

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2193-01	MGAX60	1	1.11	9.66	7.55	75.3	
L2193-02	MGAX61	2	1.13	9.68	7.73	77.2	
L2193-03	MGAX62	3	1.19	9.57	7.28	72.7	
L2193-04	MGAX62D	4	1.19	9.57	7.28	72.7	
L2193-05	MGAX62S	5	1.19	9.57	7.28	72.7	
L2193-06	MGAX63	6	1.12	9.67	6.77	66.1	
L2193-07	MGAX64	7	1.19	9.88	7.43	71.8	
L2193-08	MGAX65	8	1.17	9.82	7.92	78.0	
L2193-09	MGAX66	9	1.17	9.63	7.6	76.0	
L2193-10	MGAX67	10	1.1	9.72	8.7	88.2	
L2193-11	MGAX68	11	1.14	9.96	8.1	78.9	
L2193-12	MGAX69	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$