



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 04/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108670

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2122-06	SU-04-041020	1	1.14	9.9	9.44	94.7	
L2122-07	SU-04-041020	2	1.13	9.85	9.29	93.6	
L2122-08	SU-04-041020-1	3	1.15	9.54	9.31	97.3	
L2122-09	SU-04-041020-2	4	1.16	9.79	9.24	93.6	
L2229-03	TP-39	5	1.13	9.93	9.44	94.4	
L2229-07	TP-31	6	1.12	9.56	4.18	36.3	
L2229-11	TP-40	7	1.14	9.92	8.41	82.8	
L2229-19	TP-25	8	1.13	9.62	7.32	72.9	
L2229-23	TP-17	9	1.15	9.82	8.16	80.9	
L2229-27	TP-9	10	1.14	9.74	8.4	84.4	
L2229-31	TP-1	11	1.11	9.94	8.45	83.1	
L2229-35	TP-2	12	1.13	9.62	7.43	74.2	
L2230-03	TP-3	13	1.16	9.88	8.18	80.5	
L2230-07	TP-4	14	1.16	9.62	8.13	82.4	
L2230-11	TP-11	15	1.15	9.84	8.23	81.5	
L2230-27	TP-19	16	1.11	9.96	8.53	83.8	
L2230-35	TP-41	17	1.13	9.85	8.44	83.8	
L2230-43	TP-42	18	1.13	9.73	8.45	85.1	
L2242-01	B-1	19	1.15	9.93	8.29	81.3	
L2242-02	B-1-A	20	1.1	9.78	8.76	88.2	
L2242-03	B-1-B	21	1.13	9.78	8.76	88.2	
L2243-01	TP-34-33	22	1.14	7.44	6.67	87.8	
L2243-02	TP-34-33-3.5-4.0	23	1.13	9.55	8.46	87.1	
L2243-03	TP-34-33-3.0-3.5	24	1.14	9.96	8.8	86.8	
L2243-04	TP-34-33-2.5-3.0	25	1.13	9.8	8.61	86.3	
L2243-05	TP-34-33-COMP	26	1.18	9.91	8.89	88.3	
L2243-08	TP-31-30	27	1.12	5.77	5.24	88.6	
L2243-09	TP-31-30-3.5-4.0	28	1.17	9.96	8.86	87.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/13/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 04/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108670

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2243-10	TP-31-30-3.0-3.5	29	1.13	9.85	8.8	88.0	
L2243-11	TP-31-30-2.5-3.0	30	1.1	9.9	8.89	88.5	
L2243-12	TP-31-30-COMP	31	1.13	9.59	8.69	89.4	
L2243-15	TP-28	32	1.15	5.78	5.06	84.4	
L2243-16	TP-28-7.5-8.0	33	1.18	9.66	8.06	81.1	
L2243-17	TP-28-5.5-6.0	34	1.15	9.98	8.68	85.3	
L2243-18	TP-28-3.5-4.0	35	1.11	9.58	8.00	81.3	
L2243-19	TP-28-COMP	36	1.15	9.96	8.49	83.3	
L2244-01	TP-25-24	37	1.15	6.31	5.53	84.9	
L2244-02	TP-25-24-7.5-8.0	38	1.13	9.83	8.51	84.8	
L2244-03	TP-25-24-5.5-6.0	39	1.15	9.96	8.73	86.0	
L2244-04	TP-25-24-3.5-4.0	40	1.16	9.74	7.98	79.5	
L2244-05	TP-25-24-COMP	41	1.15	9.95	8.49	83.4	
L2244-08	TP-23-22	42	1.15	7.25	6.38	85.7	
L2244-09	TP-23-22-7.5-8.0	43	1.19	9.69	8.47	85.6	
L2244-10	TP-23-22-5.5-6.0	44	1.14	9.95	8.65	85.2	
L2244-11	TP-23-22-3.5-4.0	45	1.16	9.89	8.66	85.9	
L2244-12	TP-23-22-COMP	46	1.14	9.64	8.39	85.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$