



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/14/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108684

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2115-01	CL-01-04102020	1	1.15	9.76	9.08	92.1	
L2115-02	CL-01-04102020	2	1.15	9.91	9.05	90.2	
L2115-03	CL-01-04102020-1	3	1.11	9.58	9.08	94.1	
L2115-04	CL-01-04102020-2	4	1.19	9.84	9.22	92.8	
L2115-05	CL-01-04102020-3	5	1.19	9.52	8.96	93.3	
L2115-06	CL-01-04102020-4	6	1.11	9.92	9.26	92.5	
L2115-07	CL-01-04102020-5	7	1.12	9.57	8.69	89.6	
L2198-01	TP-43	8	1.18	7.24	6.27	84.0	
L2198-02	TP-43-3.0-3.5	9	1.18	9.72	8.52	85.9	
L2198-03	TP-43-2.5-3.0	10	1.15	9.97	8.8	86.7	
L2198-04	TP-43-2.0-2.5	11	1.16	9.94	8.77	86.7	
L2198-05	TP-43-COMP	12	1.18	9.54	8.58	88.5	
L2198-08	TP-42-41	13	1.11	6.52	5.88	88.2	
L2198-09	TP-42-41-8.0-8.5	14	1.18	9.69	7.96	79.7	
L2198-10	TP-42-41-6.0-6.5	15	1.16	9.71	8.39	84.6	
L2198-11	TP-42-41-4.0-4.5	16	1.18	9.98	8.8	86.6	
L2198-12	TP-42-41-COMP	17	1.14	9.96	8.45	82.9	
L2245-01	TP-5	18	1.15	9.9	7.26	69.8	
L2245-02	TP-13	19	1.16	9.75	4.08	34.0	
L2245-03	TP-7	20	1.16	9.58	7.96	80.8	
L2245-04	TP-15	21	1.18	9.58	6.54	63.8	
L2245-05	TP-14	22	1.16	9.59	7.14	70.9	
L2245-06	TP-26	23	1.14	9.94	8.92	88.4	
L2245-07	TP-26-EPH-2	24	1.16	9.77	8.22	82.0	
L2245-08	TP-26	25	1.14	9.92	8.92	88.6	
L2245-10	TP-23B	26	1.15	9.95	7.96	77.4	
L2245-11	TP-23B-EPH-2	27	1.12	9.66	7.62	76.1	
L2245-12	TP-23B	28	1.15	9.51	6.69	66.3	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/14/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108684

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2245-14	TP-23A	29	1.15	9.98	8.47	82.9	
L2245-15	TP-23A-EPH-2	30	1.11	9.61	8.08	82.0	
L2245-16	TP-23A	31	1.15	9.5	8.15	83.8	
L2245-18	TP-22	32	1.15	9.77	7.93	78.7	
L2245-19	TP-22-EPH-2	33	1.11	9.81	7.98	79.0	
L2245-20	TP-22	34	1.19	9.6	7.48	74.8	
L2245-22	TP-21A	35	1.19	9.96	8.2	79.9	
L2245-23	TP-21A-EPH-2	36	1.17	9.9	8.52	84.2	
L2245-24	TP-21A	37	1.15	9.65	8.09	81.6	
L2245-26	TP-21B	38	1.17	9.66	7.99	80.3	
L2245-27	TP-21B-EPH-2	39	1.15	9.88	8.32	82.1	
L2245-28	TP-21B	40	1.16	9.79	8.17	81.2	
L2245-30	TP-28	41	1.17	9.55	7.22	72.2	
L2245-31	TP-28-EPH-2	42	1.16	9.97	6.85	64.6	
L2245-32	TP-28	43	1.15	9.95	6.09	56.1	
L2245-34	TP-27	44	1.15	9.63	8.28	84.1	
L2245-35	TP-27-EPH-2	45	1.14	9.97	8.53	83.7	
L2245-36	TP-27	46	1.12	9.94	7.3	70.1	
L2245-38	TP-5-13-14	47	1.13	9.96	7.69	74.3	
L2245-39	TP-5-13-14-EPH-2	48	1.17	9.55	7.24	72.4	
L2245-41	TP-7-15	49	1.15	9.88	7.43	71.9	
L2245-42	TP-7-15-EPH-2	50	1.12	9.85	7.72	75.6	
L2255-01	TP-27-26	57	1.14	4.02	3.65	87.2	
L2255-02	TP-27-26-8.0-8.5	58	1.17	9.67	8.59	87.3	
L2255-03	TP-27-26-6.0-6.5	59	1.16	9.66	8.58	87.3	
L2255-04	TP-27-26-4.0-4.5	60	1.1	9.75	8.62	86.9	
L2255-05	TP-27-26-COMP	61	1.15	9.66	8.53	86.7	
L2258-01	1	51	1.14	9.51	7.39	74.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/14/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB108684

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2258-02	2	52	1.17	9.9	7.38	71.1	
L2258-03	3	53	1.16	9.56	7.6	76.7	
L2258-04	4	54	1.17	9.64	7.73	77.4	
L2258-05	5	55	1.17	9.95	8.83	87.2	
L2258-06	6	56	1.16	9.98	8.18	79.6	
L2263-01	400407703.0	62	1.15	9.95	1.62	5.3	sludge sample
L2264-01	400407721.0	63	1.15	8.52	3.29	29.0	sludge sample
L2265-01	400408020.0	64	1.13	9.88	1.98	9.7	sludge sample
L2266-01	400408026.0	65	1.11	9.71	1.58	5.5	sludge sample
L2267-01	400409714.2	66	1.11	9.97	1.39	3.2	sludge sample
L2268-01	400409717.0	67	1.12	9.88	1.81	7.9	sludge sample
L2269-01	400409752.0	68	1.18	9.94	2.82	18.7	sludge sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$