



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 5/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:35
In Date: 05/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 08:55
Out Date: 05/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109005

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2464-01	TRE-50006	1	1.00	2.00	2.00	100.0	OILY-DEBRIS,100% SOLIDS
L2466-01	DRUM-1-3	2	1.00	2.00	2.00	100.0	OILY-DEBRIS,100% SOLIDS
L2481-01	HSM-01-W	3	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-02	HSM-01-R	4	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-03	HSM-01-C	5	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-04	HSM-02-W	6	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-05	HSM-02-R	7	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-06	HSM-02-C	8	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-07	HSM-03-W	9	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-08	HSM-03-R	10	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-09	HSM-03-C	11	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-10	P-03-R	12	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-11	HSM-04-W	13	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-12	HSM-04-R	14	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-13	HSM-04-C	15	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-14	HSM-05-W	16	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-15	HSM-05-R	17	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-16	HSM-05-C	18	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-17	HSM-06-W	19	1.00	2.00	2.00	100.0	WOOD SAMPLE,100% SOLIDS
L2481-18	HSM-06-R	20	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2481-19	HSM-06-C	21	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2481-20	P-03-C	22	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-01	HSA-01-R	23	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-03	HSA-01-C	24	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-05	HSA-02-R	25	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-07	HSA-02-C	26	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-09	LP-01-R	27	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-11	LP-01-C	28	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: Daniel
Date: 5/5/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:35
In Date: 05/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 08:55
Out Date: 05/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109005

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2482-13	LP-02-R	29	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-15	LP-02-C	30	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-17	LP-03-R	31	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-19	LP-03-C	32	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-21	LW-01-R	33	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-23	LW-01-C	34	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-25	LW-02-R	35	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-27	LW-02-C	36	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-29	LW-03-R	37	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-31	LW-03-C	38	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-33	P-01-R	39	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-35	P-01-C	40	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2482-37	P-02-R	41	1.00	2.00	2.00	100.0	RUBBER SAMPLE ,100% SOLIDS
L2482-39	P-02-C	42	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2484-01	OILY-STONE	44	1.00	2.00	2.00	100.0	OILY-STONE SAMPLE,100% SOLIDS
L2486-01	400428026.0	43	1.12	9.97	1.7	6.6	sludge sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$