



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 05/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109024

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2114-01	PEN-DRUM	15	1.00	2.00	2.00	100.0	stone sample,100% solids
L2114-04	PEN-DRUM	16	1.00	2.00	2.00	100.0	stone sample,100% solids
L2492-13	OILY-DEBRIS-COMP	1	1.18	9.56	7.47	75.1	
L2497-01	HD-01-050620	17	1.15	9.63	9.08	93.5	
L2497-02	HD-01-050620	18	1.13	9.7	8.82	89.7	
L2497-03	HD-01-050620-1	19	1.12	9.6	9.22	95.5	
L2498-01	CL-02-050620	20	1.18	9.59	8.95	92.4	
L2498-02	CL-02-050620	21	1.13	9.77	8.91	90.0	
L2498-03	CL-02-050620-1	22	1.1	9.78	9.21	93.4	
L2498-04	CL-02-050620-2	23	1.12	9.89	9.48	95.3	
L2498-05	CL-02-050620-3	24	1.18	9.51	8.92	92.9	
L2498-06	CL-02-050620-4	25	1.13	9.72	9.04	92.1	
L2498-07	CL-02-050620-5	26	1.15	9.52	8.82	91.6	
L2508-06	BP-VPB180-GW-838-840	7	1.15	9.98	2.61	16.5	sludge sample
L2513-01	BU-03-050520	2	1.14	9.88	8.37	82.7	
L2513-02	BU-03-050520	3	1.17	9.72	8.14	81.5	
L2515-01	TP-18- (1-2)	4	1.18	9.89	9.06	90.5	
L2515-02	TP-18A	5	1.17	9.66	8.95	91.6	
L2515-03	TP-18B	6	1.17	9.92	9.09	90.5	
L2523-01	HDD	8	1.18	9.76	7.65	75.4	
L2523-02	HDD-A	9	1.16	9.87	7.69	75.0	
L2523-03	HDD-B	10	1.11	9.9	7.58	73.6	
L2525-01	DRUM-1	27	1.14	9.1	6.49	67.2	
L2525-02	DRUM-1-3	28	1.13	9.58	7.92	80.4	
L2529-01	72-12002	11	1.15	9.95	9.02	89.4	
L2529-02	72-12002	12	1.16	9.72	8.89	90.3	
L2529-03	72-11936	13	1.13	9.66	8.66	88.3	
L2529-04	72-11936	14	1.14	9.6	8.41	85.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/7/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 05/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109024

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2530-01	RO-14	29	1.16	9.97	9.17	90.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$