



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109054

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2495-01	HR-05-050720	26	1.12	9.77	8.9	89.9	
L2495-02	HR-05-050720	27	1.19	9.61	8.54	87.3	
L2495-03	HR-05-050720-1	28	1.13	9.5	8.53	88.4	
L2495-04	HR-05-050720-2	29	1.13	9.88	8.78	87.4	
L2503-01	OR-03-050720	30	1.15	9.64	9.09	93.5	
L2503-02	OR-03-050720	31	1.17	9.75	8.62	86.8	
L2503-03	OR-03-050720-1	32	1.16	9.83	9.52	96.4	
L2503-04	OR-03-050720-2	33	1.11	9.69	8.65	87.9	
L2506-01	TR-04-050720	34	1.16	9.64	9.04	92.9	
L2506-02	TR-04-050720	35	1.14	9.5	8.85	92.2	
L2506-03	TR-04-050720-1	36	1.15	9.89	9.34	93.7	
L2506-04	TR-04-050720-2	37	1.14	9.54	8.86	91.9	
L2531-01	TP-14- (1-2)	1	1.16	9.88	9.59	96.7	
L2531-02	TP-14A	2	1.16	9.58	8.78	90.5	
L2531-03	TP-14B	3	1.11	9.8	9.00	90.8	
L2531-05	TP-15- (1-2)	4	1.13	9.59	9.2	95.4	
L2531-06	TP-15A	5	1.15	9.83	9.16	92.3	
L2531-07	TP-15B	6	1.11	9.72	8.79	89.2	
L2533-01	MONOPOLE-FOUNDATION	7	1.00	2.00	2.00	100.0	stone sample,100% solids
L2534-01	356	8	1.00	2.00	2.00	100.0	stone sample,100% solids
L2535-01	COMP-02-06	10	1.16	9.89	9.14	91.4	
L2535-03	COMP-07-09	11	1.15	9.96	9.31	92.6	
L2540-01	ETGI-346	38	1.18	9.65	8.72	89.0	
L2540-02	ETGI-346	39	1.11	9.66	8.64	88.1	
L2540-03	HAR-STOCK	40	1.18	9.78	8.71	87.6	
L2540-04	HAR-STOCK	41	1.14	9.94	9.18	91.4	
L2540-06	HAR-STOCK-EPH-2	42	1.14	9.94	9.18	91.4	
L2540-07	ETGI-346	43	1.11	9.66	8.64	88.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/07/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/08/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109054

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2546-01	BP-VPB180-GW-918-920	9	1.19	9.98	3.14	22.2	sludge sample
L2553-01	VNJ-217	12	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2553-02	VNJ-217	13	1.00	2.00	2.00	100.0	CONCRETE SAMPLE,100% SOLIDS
L2556-01	NORTH-AOC14-CAR-WASH-TPH-1	14	1.12	9.71	8.04	80.6	
L2556-02	NORTH-AOC14-CAR-WASH-TPH-2	15	1.17	9.69	8.09	81.2	
L2556-03	NORTH-AOC14-CAR-WASH-TPH-3	16	1.14	9.96	8.33	81.5	
L2556-04	AOC14-CARWASH-TPH-4	17	1.18	9.54	7.36	73.9	
L2556-05	AOC14-CARWASH-TPH-5	18	1.13	9.78	7.65	75.4	
L2556-06	AOC14-CARWASH-TPH-6	19	1.11	9.69	7.62	75.9	
L2557-01	TP-12-(4-5)	20	1.13	9.96	8.16	79.6	
L2557-02	TP-12A	21	1.15	9.96	8.83	87.2	
L2557-03	TP12B	22	1.13	9.96	9.00	89.1	
L2557-05	TP-13-(1-2)	23	1.19	9.92	9.46	94.7	
L2557-06	TP-13A	24	1.17	9.54	8.8	91.2	
L2557-07	TP-13B	25	1.12	9.53	8.72	90.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$