

LB109109

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : du

Instrument ID : Konelab

5/12/2020 12:58

Test: ~~Total~~ ^{Relative} CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	101.203	0.0	0.062	
ICB1	-0.268	0.0	0.003	
CCV1	255.100	0.0	0.150	
CCB1	-0.261	0.0	0.003	
PB128817BL	-0.612	0.0	0.003	
L2495-05	-0.622	0.0	0.003	
L2495-05DUP	-0.688	0.0	0.003	
L2502-06	-0.675	0.0	0.003	
L2503-05	-0.152	0.0	0.003	
L2557-04	-0.371	0.0	0.003	
L2557-08	-0.225	0.0	0.003	
L2571-01	-0.496	0.0	0.003	
L2571-02	-0.324	0.0	0.003	
L2571-03	-0.748	0.0	0.003	
CCV2	273.468	0.0	0.161	
CCB2	-0.643	0.0	0.003	
L2571-04	-0.814	0.0	0.003	
L2571-05	-0.963	0.0	0.003	
PB128831BL	-0.424	0.0	0.003	
L2576-04	-0.875	0.0	0.003	
L2576-04DUP	-0.994	0.0	0.003	
L2531-04	-0.888	0.0	0.003	
L2531-08	-0.566	0.0	0.003	
L2554-02	-0.800	0.0	0.003	
L2554-04	-0.240	0.0	0.003	
CCV3	254.170	0.0	0.149	
CCB3	-0.600	0.0	0.003	

N	27
Mean	32.248
SD	84.7250
CV%	262.73

Aquakem v. 7.2AQ1

Results from time period:

Tue May 12 10:05:02 2020

Tue May 12 12:56:50 2020

174

Sample Id	Sam/Ctr/c/	Test short nam	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	REACTIVE CN	P	-1.0828	µg/l	5/12/2020 11:26:00	
5.OPPBCN	A	REACTIVE CN	P	4.7695	µg/l	5/12/2020 11:26:01	
10PPBCN	A	REACTIVE CN	P	10.6678	µg/l	5/12/2020 11:26:02	
50PPBCN	A	REACTIVE CN	P	50.7759	µg/l	5/12/2020 11:26:03	
100PPBCN	A	REACTIVE CN	P	100.9102	µg/l	5/12/2020 11:26:04	
250PPBCN	A	REACTIVE CN	P	248.4601	µg/l	5/12/2020 11:26:05	
500PPBCN	A	REACTIVE CN	P	500.4992	µg/l	5/12/2020 11:26:06	
ICV1	S	REACTIVE CN	P	101.2026	µg/l	5/12/2020 12:16:21	
ICB1	S	REACTIVE CN	P	-0.2679	µg/l	5/12/2020 12:16:24	
CCV1	S	REACTIVE CN	P	255.1001	µg/l	5/12/2020 12:16:26	
CCB1	S	REACTIVE CN	P	-0.2613	µg/l	5/12/2020 12:16:27	
PB128817BL	S	REACTIVE CN	P	-0.6122	µg/l	5/12/2020 12:16:30	
L2495-05	S	REACTIVE CN	P	-0.6225	µg/l	5/12/2020 12:16:31	
L2495-05DUP	S	REACTIVE CN	P	-0.6885	µg/l	5/12/2020 12:23:53	
L2502-06	S	REACTIVE CN	P	-0.6754	µg/l	5/12/2020 12:23:54	
L2503-05	S	REACTIVE CN	P	-0.1516	µg/l	5/12/2020 12:23:55	
L2557-04	S	REACTIVE CN	P	-0.3713	µg/l	5/12/2020 12:23:56	
L2557-08	S	REACTIVE CN	P	-0.2252	µg/l	5/12/2020 12:23:57	
L2571-01	S	REACTIVE CN	P	-0.4963	µg/l	5/12/2020 12:23:58	
L2571-02	S	REACTIVE CN	P	-0.3237	µg/l	5/12/2020 12:23:59	
L2571-03	S	REACTIVE CN	P	-0.7484	µg/l	5/12/2020 12:24:00	
CCV2	S	REACTIVE CN	P	273.4679	µg/l	5/12/2020 12:24:01	
CCB2	S	REACTIVE CN	P	-0.6434	µg/l	5/12/2020 12:28:40	
L2571-04	S	REACTIVE CN	P	-0.8135	µg/l	5/12/2020 12:28:41	
L2571-05	S	REACTIVE CN	P	-0.9632	µg/l	5/12/2020 12:28:42	
PB128831BL	S	REACTIVE CN	P	-0.4241	µg/l	5/12/2020 12:28:44	
L2576-04	S	REACTIVE CN	P	-0.8747	µg/l	5/12/2020 12:34:44	
L2576-04DUP	S	REACTIVE CN	P	-0.9944	µg/l	5/12/2020 12:34:46	
L2531-04	S	REACTIVE CN	P	-0.8885	µg/l	5/12/2020 12:34:47	
L2531-08	S	REACTIVE CN	P	-0.5659	µg/l	5/12/2020 12:34:48	
L2554-02	S	REACTIVE CN	P	-0.7997	µg/l	5/12/2020 12:39:59	
L2554-04	S	REACTIVE CN	P	-0.2398	µg/l	5/12/2020 12:40:00	
CCV3	S	REACTIVE CN	P	254.17	µg/l	5/12/2020 12:56:47	
CCB3	S	REACTIVE CN	P	-0.5998	µg/l	5/12/2020 12:56:50	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : De

Instrument ID : Konelab

5/12/2020 11:31


Test ^{Leucine} ~~Total~~ CN

Accepted 5/12/2020 11:31

Factor 1740
Bias 0.003

Coeff. of det. 0.999972

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	-1.0828	0.0000	
2	5.0PPBCN	0.006	4.7695	5.0000	
3	10PPBCN	0.009	10.6678	10.0000	
4	50PPBCN	0.033	50.7759	50.0000	
5	100PPBCN	0.061	100.9102	100.0000	
6	250PPBCN	0.146	248.4601	250.0000	
7	500PPBCN	0.291	500.4992	500.0000	