



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/14/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 05/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 05/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109143

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2501-01	EO-01-051320	1	1.13	9.66	8.99	92.1	
L2501-02	EO-01-051320	2	1.17	9.59	9.03	93.3	
L2501-03	EO-01-051320-1	3	1.13	9.72	9.03	92.0	
L2501-04	EO-01-051320-2	4	1.11	9.66	8.71	88.9	
L2501-05	EO-01-051320-3	5	1.15	9.55	8.86	91.8	
L2501-06	EO-01-051320-4	6	1.18	9.79	9.52	96.9	
L2616-01	NOT-APPLICABLE	7	1.00	2.00	2.00	100.0	
L2635-01	ROLL-OFF-P72-11933	8	1.12	9.61	8.76	90.0	
L2636-01	AST-1	9	1.18	9.88	8.55	84.7	
L2636-02	AST-2	10	1.19	9.98	8.7	85.4	
L2636-03	AST-3	11	1.18	9.94	8.57	84.4	
L2638-02	SUP-WZ2-G-1	12	1.11	9.5	9.19	96.3	
L2638-03	SUP-WZ2-G-2	13	1.12	9.77	8.7	87.6	
L2638-04	SUP-WZ2-G-3	14	1.18	9.66	9.00	92.2	
L2638-05	SUP-WZ2-C-1	15	1.19	9.7	8.82	89.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$