



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/28/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 16:30
In Date: 05/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 05/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109330

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2505-06	SU-01-052720	40	1.12	9.87	9.28	93.3	
L2505-07	SU-01-052720	41	1.14	9.75	9.00	91.3	
L2505-08	SU-01-052720-1	42	1.11	9.96	8.68	85.5	
L2505-09	SU-01-052720-2	43	1.18	9.84	9.57	96.9	
L2767-01	TRE-1005	1	1.00	2.00	2.00	100.0	OILY DEBRIS,100% SOLIDS
L2773-01	NM10-COMP-2020-0-6	2	1.15	9.77	8.27	82.6	
L2773-02	NM8-COMP-2020-0-6	3	1.14	9.56	8.18	83.6	
L2773-03	L2773-02MS	4	1.14	9.56	8.18	83.6	
L2773-04	L2773-02MSD	5	1.14	9.56	8.18	83.6	
L2773-05	NM16-COMP-2020-0-6	6	1.18	9.8	8.38	83.5	
L2773-06	NMDUP04-COMP-2020-0-6	7	1.14	9.94	8.67	85.6	
L2773-07	NM74-COMP-2020-0-6	8	1.18	9.64	8.24	83.5	
L2773-08	NM59-COMP-2020-0-6	9	1.16	9.97	7.6	73.1	
L2773-09	NM35-COMP-2020-0-6	10	1.17	9.93	7.51	72.4	
L2773-10	NM46-COMP-2020-0-6	11	1.12	9.84	7.39	71.9	
L2773-11	NM47-COMP-2020-0-6	12	1.19	9.67	7.66	76.3	
L2773-12	NM48-COMP-2020-0-6	13	1.16	9.55	8.32	85.3	
L2773-13	NM4-COMP-2020-0-6	14	1.13	9.52	7.88	80.5	
L2773-14	NM25-COMP-2020-0-6	15	1.16	9.9	8.27	81.4	
L2773-15	NM26-COMP-2020-0-6	16	1.17	9.54	7.86	79.9	
L2773-16	NM61-COMP-2020-0-6	17	1.18	9.62	8.24	83.6	
L2773-17	NM33-COMP-2020-0-6	18	1.15	9.96	8.99	89.0	
L2773-19	NM58-COMP-2020-0-6	19	1.17	9.68	7.37	72.9	
L2773-20	NM73-COMP-2020-0-6	20	1.13	9.94	7.05	67.2	
L2773-21	NM91-COMP-2020-0-6	21	1.17	9.7	8.44	85.2	
L2773-22	NM92-COMP-2020-0-6	22	1.13	9.92	8.74	86.6	
L2773-23	NM89-COMP-2020-0-6	23	1.15	9.57	8.12	82.8	
L2773-24	NM90-COMP-2020-0-6	24	1.13	9.86	8.00	78.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/28/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 16:30
In Date: 05/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 05/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109330

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2773-25	NM113-COMP-2020-0-6	25	1.1	9.83	8.27	82.1	
L2773-26	NM3-COMP-2020-0-6	26	1.16	9.93	8.9	88.3	
L2775-01	PAINTED-CURB	44	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2775-02	UNPAINTED-CURB	45	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2776-01	351	27	1.14	9.62	6.58	64.2	
L2776-02	351	28	1.19	9.8	6.6	62.8	
L2776-03	225	29	1.13	9.84	8.4	83.5	
L2776-04	225	30	1.15	9.98	8.73	85.8	
L2776-05	339	31	1.12	9.52	5.52	52.4	
L2776-06	339	32	1.1	9.76	6.28	59.8	
L2783-01	LOT-11	33	1.17	9.61	8.95	92.2	
L2784-01	B-1	34	1.14	9.7	8.59	87.0	
L2784-02	B-2	35	1.18	9.96	8.83	87.1	
L2784-03	B-3	36	1.12	9.93	8.71	86.2	
L2784-04	B-4	37	1.19	9.64	8.84	90.5	
L2784-05	B-5	38	1.18	9.98	9.32	92.5	
L2784-06	B-6	39	1.12	9.52	8.8	91.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$