



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/29/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2503-11	OR-03-052820	30	1.14	9.58	9.02	93.4	
L2503-12	OR-03-052820	31	1.15	9.96	9.45	94.2	
L2503-13	OR-03-052820-1	32	1.17	9.98	9.7	96.8	
L2503-14	OR-03-052820-2	33	1.13	9.94	9.12	90.7	
L2540-08	HAR-STOCK	1	1.14	9.94	9.18	91.4	
L2790-01	SB-1-S1	15	1.13	9.64	8.58	87.5	
L2790-02	SB-1-S2	16	1.15	9.97	9.00	89.0	
L2790-03	SB-1-S	17	1.12	9.65	8.48	86.3	
L2790-05	SB-1-D1	18	1.13	9.88	8.84	88.1	
L2790-06	SB-1-D2	19	1.15	9.96	9.24	91.8	
L2790-07	SB-1-D	20	1.12	9.72	8.7	88.1	
L2790-09	SB-3-S1	21	1.17	9.68	9.17	94.0	
L2790-10	SB-3-S2	22	1.1	9.58	9.09	94.2	
L2790-11	SB-3-S	23	1.15	9.68	9.18	94.1	
L2790-15	SB-7-S1	24	1.16	9.93	9.16	91.2	
L2790-16	SB-7-S2	25	1.13	9.77	9.1	92.2	
L2790-17	SB-7-S	26	1.19	9.61	9.02	93.0	
L2792-01	NM32-COMP-2020-0-6	2	1.12	9.83	8.02	79.2	
L2792-02	NM28-COMP-2020-0-6	3	1.17	9.58	8.41	86.1	
L2792-03	L2792-02MS	4	1.17	9.58	8.41	86.1	
L2792-04	L2792-02MSD	5	1.17	9.58	8.41	86.1	
L2792-05	NM15-COMP-2020-0-6	6	1.13	9.8	8.32	82.9	
L2792-06	NM14-COMP-2020-0-6	7	1.18	9.56	8.23	84.1	
L2792-07	NM107-COMP-2020-0-6	8	1.11	9.95	8.49	83.5	
L2792-08	NM114-COMP-2020-0-6	9	1.17	9.96	8.85	87.4	
L2792-09	NM84-COMP-2020-0-6	10	1.17	9.82	7.16	69.2	
L2792-10	NM96-COMP-2020-0-6	11	1.14	9.51	7.54	76.5	
L2792-11	NM85-COMP-2020-0-6	12	1.18	9.64	8.18	82.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/29/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109343

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2792-12	NM75-COMP-2020-0-6	13	1.17	9.69	8.06	80.9	
L2792-13	NM97-COMP-2020-0-6	14	1.16	9.83	7.95	78.3	
L2795-01	217	27	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2795-03	209-213	28	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2795-04	PALLET	29	1.00	2.00	2.00	100.0	WOOD CHIPS SAMPLE,100% SOLIDS
L2796-01	CON-WEST-A	34	1.00	2.00	2.00	100.0	stone sample,100% solids
L2796-02	CON-WEST-B	35	1.00	2.00	2.00	100.0	stone sample,100% solids
L2796-03	CON-EAST-A	36	1.00	2.00	2.00	100.0	stone sample,100% solids
L2796-04	CON-EAST-B	37	1.00	2.00	2.00	100.0	stone sample,100% solids
L2796-05	CON-EAST-C	38	1.00	2.00	2.00	100.0	stone sample,100% solids
L2797-01	TEST-HOLE-5-COMP	39	1.13	9.62	8.18	83.0	
L2797-04	TEST-HOLE-5-EPH-1.5-2.0	40	1.14	9.9	8.46	83.6	
L2797-05	TEST-HOLE-5-VOC-1.5-2.0	41	1.18	9.94	8.49	83.4	
L2797-06	TEST-HOLE-3-COMP	42	1.13	9.55	8.64	89.2	
L2797-09	TEST-HOLE-3-EPH-2.0-2.5	43	1.11	9.79	9.12	92.3	
L2797-10	TEST-HOLE-3-VOC-2.0-2.5	44	1.12	9.82	9.15	92.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$