



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/4/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 06/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 06/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109417

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2877-01	ETX13	1	1.13	9.58	8.5	87.2	
L2877-02	ETX14	2	1.17	9.62	8.68	88.9	
L2877-03	ETX14MS	3	1.17	9.62	8.68	88.9	
L2877-04	ETX14MSD	4	1.17	9.62	8.68	88.9	
L2877-05	ETX15	5	1.13	9.82	8.95	90.0	
L2877-06	ETX16	6	1.11	9.92	9.06	90.2	
L2877-07	ETX17	7	1.17	9.58	8.54	87.6	
L2877-08	ETX18	8	1.1	9.74	8.88	90.0	
L2877-09	ETX19	9	1.15	9.55	8.7	89.9	
L2877-10	ETX20	10	1.16	9.86	9.08	91.0	
L2877-11	ETX21	11	1.18	9.68	8.87	90.5	
L2877-12	ETX22	12	1.12	9.72	9.00	91.6	
L2877-13	ETX23	13	1.15	9.6	8.84	91.0	
L2877-14	ETX24	14	1.14	9.9	8.81	87.6	
L2877-15	ETX25	15	1.15	9.64	8.73	89.3	
L2877-16	ETX26	16	1.14	9.56	8.73	90.1	
L2877-17	ETX27	17	1.13	9.66	8.56	87.1	
L2877-18	ETX29	18	1.13	9.66	8.24	83.4	
L2877-19	ETX30	19	1.15	9.95	9.83	98.6	
L2877-20	ETX31	20	1.19	9.77	7.87	77.9	
L2877-21	ETX32	21	1.12	9.97	8.67	85.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$