



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/5/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 16:35
In Date: 06/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 06/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109440

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2837-01	SP	1	1.12	9.84	9.43	95.3	
L2886-01	SLAB-4	11	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2886-02	GEN	12	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2886-04	SLAB-3	13	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2886-05	SLAB-2	14	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2886-06	SLAB-1	15	1.00	2.00	2.00	100.0	STONE SAMPLE,100% solids
L2891-01	PCB-060120-SR-03	2	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2892-01	PCB-060120-SR-02	3	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2893-01	SP-2A	4	1.19	9.98	9.51	94.7	
L2918-01	1-19029-TC-1	5	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L2918-02	2-19029-TC-1	6	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L2918-03	3-19029-TC-1	7	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L2918-04	4-19029-TC-1	8	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L2918-05	5-19029-TC-1	9	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L2918-06	6-19029-TC-1	10	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$