



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109450

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2812-01	OK-01-060520	60	1.14	9.96	9.34	93.0	
L2812-02	OK-01-060520	61	1.14	9.78	8.85	89.2	
L2812-03	OK-01-060520-1	62	1.16	9.97	9.26	91.9	
L2812-04	OK-01-060520-2	63	1.13	9.71	9.27	94.9	
L2815-01	OR-02-060520	64	1.18	9.9	8.96	89.2	
L2815-02	OR-02-060520	65	1.11	9.64	8.82	90.4	
L2815-03	OR-02-060520-1	66	1.17	9.7	8.61	87.2	
L2815-04	OR-02-060520-2	67	1.13	9.55	8.51	87.6	
L2887-01	CHRT-25606	1	1.15	9.98	5.45	48.7	SLUDGE SAMPLE
L2887-03	ETGI-335	2	1.16	9.81	9.45	95.8	
L2890-01	CHRT-26160	3	1.19	9.75	9.06	91.9	
L2890-03	72-11979-72-11969	4	1.13	9.63	8.82	90.5	
L2890-05	CHRT-24245	5	1.14	9.67	6.75	65.8	
L2890-07	CHRT-26572-CHRT-25062-R BR-251686	6	1.12	9.62	7.39	73.8	
L2922-05	SVOC-GPC-BLANK	26	1.00	2.00	2.00	100.0	
L2922-06	PEST-GPC-BLANK	27	1.00	2.00	2.00	100.0	
L2922-07	PEST-GPC-BLANK-SPIKE	28	1.00	2.00	2.00	100.0	
L2922-08	PCB-GPC-BLANK	29	1.00	2.00	2.00	100.0	
L2922-09	PCB-GPC-BLANK-SPIKE	30	1.00	2.00	2.00	100.0	
L2922-10	SVOC-GCP2-BLANK	31	1.00	2.00	2.00	100.0	
L2922-11	PEST-GCP2-BLANK	32	1.00	2.00	2.00	100.0	
L2922-12	PEST-GCP2-BLANK-SPIKE	33	1.00	2.00	2.00	100.0	
L2922-13	PCB-GCP2-BLANK	34	1.00	2.00	2.00	100.0	
L2922-14	PCB-GCP2-BLANK-SPIKE	35	1.00	2.00	2.00	100.0	
L2923-01	31008	90	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2923-02	31048	91	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2923-03	COMP-2	92	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2924-01	SP	93	1.12	9.77	8.13	81.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109450

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2924-02	SP-EPH-2	94	1.12	9.77	8.13	81.0	
L2930-01	S39-HSVE2-IDW01	36	1.15	9.95	8.15	79.5	
L2930-02	S39-HSVE3-IDW01	37	1.13	9.98	8.46	82.8	
L2932-01	PCB-042720-SR-01	7	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-02	PCB-042720-SR-02	8	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-03	PCB-042720-SR-03	9	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-04	PCB-042720-SR-04	10	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-05	PCB-042720-SR-05	11	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-06	PCB-042720-SR-06	12	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-07	PCB-042720-SR-07	13	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-08	PCB-042720-SR-08	14	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-09	PCB-042720-SR-09	15	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-10	PCB-042720-SR-10	16	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-11	PCB-042720-SR-11	17	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-12	PCB-042720-SR-12	18	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-13	PCB-042720-SR-13	19	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-14	PCB-042720-SR-14	20	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-15	PCB-042720-SR-15	21	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-16	PCB-042720-SR-16	22	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-17	PCB-042720-SR-17	23	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-18	PCB-042720-SR-18	24	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2932-19	PCB-042720-SR-19	25	1.00	2.00	2.00	100.0	CAULKING SAMPLE,100% SOLIDS
L2936-01	WASTE	38	1.14	9.71	8.8	89.4	
L2936-02	VOC	39	1.13	9.57	8.72	89.9	
L2936-03	1	40	1.14	9.95	9.28	92.4	
L2936-04	2	41	1.16	9.93	9.06	90.1	
L2936-05	3	42	1.14	9.84	8.86	88.7	
L2936-06	4	43	1.17	9.67	8.89	90.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109450

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2936-07	5	44	1.17	9.97	9.4	93.5	
L2937-01	326FB-001-0004	68	1.17	9.77	8.16	81.3	
L2937-02	326FB-002-0004	69	1.13	9.77	8.66	87.2	
L2937-03	326FB-003-0004	70	1.16	9.8	8.5	85.0	
L2937-04	L2937-03MS	71	1.16	9.8	8.5	85.0	
L2937-05	L2937-03MSD	72	1.16	9.8	8.5	85.0	
L2937-06	326FB-004-0004	73	1.14	9.9	6.16	57.3	
L2937-07	2020-06-01-FD01	74	1.1	9.74	8.00	79.9	
L2937-08	326FB-005-0004	75	1.14	9.96	6.88	65.1	
L2937-09	326FB-006-0004	76	1.15	9.98	8.73	85.8	
L2937-10	326FB-007-0004	77	1.15	9.93	8.36	82.1	
L2937-11	326FB-008-0004	78	1.1	9.85	9.00	90.3	
L2937-12	326FB-009-0004	79	1.12	9.91	8.28	81.5	
L2937-13	326FB-013-0004	80	1.13	9.96	9.49	94.7	
L2937-14	200FB-001-0004	81	1.18	9.92	8.56	84.4	
L2937-15	200FB-002-0004	82	1.13	9.9	8.85	88.0	
L2937-16	200FB-003-0004	83	1.14	9.79	7.63	75.0	
L2937-17	326FB-012-0003	84	1.14	9.96	8.33	81.5	
L2937-18	326FB-014-000105	85	1.11	9.83	8.38	83.4	
L2937-19	W00SPK-001-0004	86	1.18	9.58	7.06	70.0	
L2937-20	W00SPK-002-0004	87	1.14	9.72	8.98	91.4	
L2937-21	W00SPK-003-0004	88	1.13	9.89	8.69	86.3	
L2937-22	W00SPK-004-0004	89	1.12	9.66	8.36	84.8	
L2938-01	W00SPK-005-0004	45	1.15	9.93	9.41	94.1	
L2938-02	W00SPK-006-0004	46	1.11	9.63	8.19	83.1	
L2938-03	2020-06-02-FD03	47	1.16	9.7	8.2	82.4	
L2938-04	W00SPK-007-0004	48	1.19	9.92	9.31	93.0	
L2938-05	L2938-04MS	49	1.19	9.92	9.31	93.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/06/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109450

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2938-06	L2938-04MSD	50	1.19	9.92	9.31	93.0	
L2938-07	W00SPK-008-0004	51	1.1	9.6	8.36	85.4	
L2938-08	W00SPK-009-0004	52	1.16	9.77	7.91	78.4	
L2938-09	W00SPK-010-0004	53	1.14	9.82	9.55	96.9	
L2938-10	W00SPK-011-0004	54	1.1	9.92	7.18	68.9	
L2938-11	W00SPK-012-0004	55	1.15	9.94	9.36	93.4	
L2938-12	W00SPK-013-0004	56	1.11	9.71	8.81	89.5	
L2938-13	W00SPK-014-0004	57	1.16	9.93	9.01	89.5	
L2938-14	W00SPK-015-0004	58	1.13	9.53	7.1	71.1	
L2939-01	PHC-326FB-010	95	1.18	9.83	7.77	76.2	
L2939-02	2020-06-02-FD02	96	1.18	9.6	7.66	77.0	
L2939-03	PHC-326FB-011	97	1.1	9.81	7.84	77.4	
L2941-01	RBF-TS-001	59	1.19	9.62	7.74	77.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$