



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 16:50
In Date: 06/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 06/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109520

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2810-01	CL-02-060920	10	1.13	9.73	9.28	94.8	
L2810-02	CL-02-060920	11	1.18	9.56	9.15	95.1	
L2810-03	CL-02-060920-1	12	1.15	9.91	9.44	94.6	
L2810-04	CL-02-060920-2	13	1.16	9.82	9.07	91.3	
L2810-05	CL-02-060920-3	14	1.12	9.82	9.27	93.7	
L2810-06	CL-02-060920-4	15	1.16	9.83	9.33	94.2	
L2810-07	CL-02-060920-5	16	1.1	9.55	9.18	95.6	
L2923-13	23850	1	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L2954-01	239-355	2	1.14	9.59	8.81	90.8	
L2954-02	239-355	3	1.13	9.68	8.96	91.6	
L2954-04	239-355-EPH-2	4	1.13	9.68	8.96	91.6	
L2961-01	FM-CON-1	5	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L2961-02	FM-CON-2	6	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L2961-03	FM-CON-3	7	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L2961-04	FM-CON-4	8	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L2961-05	FM-CON-5	9	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$