



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/11/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 13:56
In Date: 06/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:42
Out Date: 06/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109543

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L2966-01	ETX99	1	1.17	9.87	9.06	90.7	
L2966-02	ETXA0	2	1.18	9.62	8.94	91.9	
L2966-03	ETXA1	3	1.17	9.53	8.78	91.0	
L2966-04	ETXA2	4	1.12	9.63	8.8	90.2	
L2966-05	ETXA3	5	1.16	9.55	8.49	87.4	
L2966-06	ETXA4	6	1.16	9.67	8.71	88.7	
L2966-07	ETXA4MS	7	1.16	9.67	8.71	88.7	
L2966-08	ETXA4MSD	8	1.16	9.67	8.71	88.7	
L2966-09	ETXA5	9	1.15	9.8	8.66	86.8	
L2966-10	ETXA6	10	1.14	9.53	7.98	81.5	
L2966-11	ETXA7	11	1.12	9.78	8.12	80.8	
L2966-12	ETXA8	12	1.16	9.98	8.88	87.5	
L2966-13	ETXA9	13	1.14	9.97	9.00	89.0	
L2966-14	ETXB0	14	1.13	9.65	7.89	79.3	
L2966-15	ETXB1	15	1.14	9.57	8.72	89.9	
L2966-16	ETXB2	16	1.17	9.81	9.16	92.5	
L2966-17	ETXB3	17	1.15	9.74	8.81	89.2	
L2966-18	ETXB4	18	1.15	9.76	8.62	86.8	
L2966-19	ETXB5	19	1.15	9.68	8.59	87.2	
L2966-20	ETXB6	20	1.15	9.75	8.65	87.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$