



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/15/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 06/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109572

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2995-01	239-WC-S-10-0-2	1	1.15	9.94	8.29	81.2	
L2995-02	239-WC-S-10-2-4	2	1.18	9.66	8.41	85.3	
L2995-03	239-WC-S-8-0-2	3	1.15	9.62	8.45	86.2	
L2995-04	239-WC-S-8-2-4	4	1.19	9.95	8.75	86.3	
L2995-05	239-WC-S-7-0-2	5	1.12	9.68	8.94	91.4	
L2995-06	239-WC-S-7-2-4	6	1.13	9.94	8.98	89.1	
L2995-07	239-WC-S-107-0-2	7	1.12	9.58	8.79	90.7	
L2995-08	239-WC-S-9-0-2	8	1.13	9.6	8.94	92.2	
L2995-09	239-WC-S-9-2-4	9	1.17	9.74	8.32	83.4	
L2995-10	239-WC-S-11-0-2	10	1.17	9.88	8.71	86.6	
L2995-11	239-WC-S-11-2-4	11	1.14	9.96	8.97	88.8	
L2995-12	239-WC-S-11-4-6	12	1.1	9.91	8.67	85.9	
L2995-13	239-WC-S-6-0-2	13	1.13	9.65	8.53	86.9	
L2995-14	239-WC-S-6-2-4	14	1.13	9.88	8.77	87.3	
L2995-15	239-WC-S-6-4-6	15	1.14	9.96	8.66	85.3	
L2995-16	239-WC-S-1-0-2	16	1.17	9.77	8.38	83.8	
L2995-17	239-WC-S-1-2-4	17	1.1	9.59	8.49	87.0	
L2995-18	239-WC-S-1-4-6	18	1.16	9.97	8.65	85.0	
L2995-19	239-WC-S-3-0-2	19	1.16	9.71	8.79	89.2	
L2995-20	239-WC-S-3-2-4	20	1.13	9.56	8.04	82.0	
L2995-21	239-WC-S-3-4-6	21	1.13	9.54	8.34	85.7	
L2995-22	239-WC-S-15-0-1.5	22	1.17	9.98	9.2	91.1	
L2995-23	239-WC-S-15-1.5-3	23	1.17	9.67	8.59	87.3	
L2995-24	239-WC-S-16-0-1.5	24	1.19	9.77	8.9	89.9	
L2995-25	239-WC-S-16-1.5-3	25	1.14	9.93	8.98	89.2	
L2996-01	DRUM-COMP	36	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L2997-01	294	37	1.13	9.95	9.00	89.2	
L2997-03	4352	38	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/15/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 06/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109572

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2998-01	BS-DRUMS	39	1.14	9.93	9.44	94.4	
L2999-03	TRE-1013	40	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L2999-05	30662	41	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L3000-05	SVOC-GPC-BLANK	26	1.00	2.00	2.00	100.0	
L3000-06	PEST-GPC-BLANK	27	1.00	2.00	2.00	100.0	
L3000-07	PEST-GPC-BLANK-SPIKE	28	1.00	2.00	2.00	100.0	
L3000-08	PCB-GPC-BLANK	29	1.00	2.00	2.00	100.0	
L3000-09	PCB-GPC-BLANK-SPIKE	30	1.00	2.00	2.00	100.0	
L3000-10	SVOC-GPC2-BLANK	31	1.00	2.00	2.00	100.0	
L3000-11	PEST-GPC2-BLANK	32	1.00	2.00	2.00	100.0	
L3000-12	PEST-GPC2-BLANK-SPIKE	33	1.00	2.00	2.00	100.0	
L3000-13	PCB-GPC2-BLANK	34	1.00	2.00	2.00	100.0	
L3000-14	PCB-GPC2-BLANK-SPIKE	35	1.00	2.00	2.00	100.0	
L3001-01	PAD	42	1.12	9.84	8.98	90.1	
L3010-01	239-WC-S-2-0-2	43	1.13	9.92	8.65	85.6	
L3010-02	239-WC-S-2-2-4	44	1.16	9.85	8.66	86.3	
L3010-03	239-WC-S-2-4-6	45	1.15	9.77	8.53	85.6	
L3010-04	239-WC-S-102-2-4	46	1.18	9.72	8.66	87.6	
L3010-05	239-WC-S-17-0-1.5	47	1.18	9.59	8.71	89.5	
L3010-06	239-WC-S-17-1.5-3	48	1.13	9.61	8.65	88.7	
L3010-07	239-WC-S-18-0-1.5	49	1.14	9.97	8.68	85.4	
L3010-08	239-WC-S-18-1.5-3	50	1.14	9.92	8.87	88.0	
L3010-09	239-WC-S-5-0-2	51	1.18	9.62	8.05	81.4	
L3010-10	239-WC-S-5-2-4	52	1.18	9.94	8.71	86.0	
L3010-11	239-WC-S-5-4-6	53	1.18	9.59	8.37	85.5	
L3010-12	239-WC-S-19-0-1.5	54	1.13	9.86	8.98	89.9	
L3010-13	239-WC-S-19-1.5-3	55	1.13	9.77	8.63	86.8	
L3010-14	239-WC-S-4-0-2.5	56	1.17	9.68	8.22	82.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/15/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 06/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109572

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3010-15	239-WC-S-4-2.5-5	57	1.13	9.97	8.33	81.4	
L3010-16	239-WC-S-12-0-2.5	58	1.12	9.72	8.7	88.1	
L3010-17	239-WC-S-12-2.5-5	59	1.14	9.97	8.76	86.3	
L3010-18	239-WC-S-13-0-2.5	60	1.14	9.97	8.94	88.3	
L3010-19	239-WC-S-13-2.5-5	61	1.14	9.96	8.74	86.2	
L3010-20	239-WC-S-14-0-1.5	62	1.14	9.79	8.56	85.8	
L3010-21	239-WC-S-14-1.5-3	63	1.15	9.98	9.02	89.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$