



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/22/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 06/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 06/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109680

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3082-01	SB-1	1	1.13	9.55	8.86	91.8	
L3082-02	SB-2	2	1.14	9.91	9.18	91.7	
L3082-03	SB-3	3	1.14	9.82	8.87	89.1	
L3082-04	SB-4	4	1.15	9.84	8.88	89.0	
L3082-05	SB-5	5	1.12	9.78	8.85	89.3	
L3082-06	SB-7	6	1.14	9.69	9.14	93.6	
L3082-07	SB-8	7	1.13	9.94	9.31	92.8	
L3082-08	SB-9	8	1.16	9.85	9.33	94.0	
L3082-09	SB-10	9	1.11	9.92	8.55	84.4	
L3082-10	SB-11	10	1.13	9.73	8.73	88.4	
L3082-11	SB-12	11	1.15	9.54	8.52	87.8	
L3082-12	SB-13	12	1.16	9.89	9.06	90.5	
L3082-13	SB-14	13	1.12	9.9	8.82	87.7	
L3082-14	SB-15	14	1.13	9.76	8.2	81.9	
L3082-15	SB-16	15	1.15	9.98	9.15	90.6	
L3082-16	SB-17	16	1.14	9.59	8.48	86.9	
L3082-17	SB-18	17	1.13	9.97	9.31	92.5	
L3082-18	SB-19	18	1.12	9.93	9.06	90.1	
L3082-19	SB-20	19	1.13	9.98	9.37	93.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$