



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 06/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109785

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L2815-16	OR-03-062920	46	1.13	9.77	9.08	92.0	
L2815-17	OR-03-062920	47	1.13	9.65	8.93	91.5	
L2815-18	OR-03-062920-1	48	1.19	9.98	9.18	90.9	
L2815-19	OR-03-062920-2	49	1.11	9.77	9.03	91.5	
L3148-01	EPEC-POLYMERS-01	1	1.12	9.75	5.67	52.7	
L3148-02	EPEC-POLYMERS-02	2	1.13	9.97	8.7	85.6	
L3148-03	EPEC-POLYMERS-03	3	1.18	9.98	4.19	34.2	
L3148-04	EPEC-POLYMERS-04	4	1.13	9.83	6.14	57.6	
L3148-05	EPEC-POLYMERS-05	5	1.16	9.98	5.82	52.8	
L3151-01	BP-MW180-SOIL-580-582	6	1.14	9.94	7.58	73.2	
L3151-02	BP-MW180-584-588	7	1.12	9.97	8.84	87.2	
L3152-02	40680	8	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3152-03	DRUM-COMP	9	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3153-01	SP1-A	10	1.12	9.73	8.61	87.0	
L3153-03	SP1-A-EPH-2	11	1.12	9.73	8.61	87.0	
L3153-04	SP1-B	12	1.16	9.88	8.54	84.6	
L3153-06	SP1-B-EPH-2	13	1.16	9.88	8.54	84.6	
L3153-07	SP1-C	14	1.17	9.55	8.22	84.1	
L3153-09	SP1-C-EPH-2	15	1.17	9.55	8.22	84.1	
L3153-10	SP2-A	16	1.15	9.95	8.64	85.1	
L3153-12	SP2-A-EPH-2	17	1.15	9.95	8.64	85.1	
L3153-13	SP2-B	18	1.17	9.73	8.5	85.6	
L3153-15	SP2-B-EPH-2	19	1.17	9.73	8.5	85.6	
L3153-16	SP2-C	20	1.16	9.79	8.43	84.2	
L3153-18	SP2-C-EPH-2	21	1.16	9.79	8.43	84.2	
L3153-19	SP2-D	22	1.12	9.9	8.51	84.2	
L3153-21	SP2-D-EPH-2	23	1.12	9.9	8.51	84.2	
L3153-22	SP2-E	24	1.19	9.75	8.59	86.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/30/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 06/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 06/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109785

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3153-24	SP2-E-EPH-2	25	1.19	9.75	8.59	86.4	
L3153-25	SP2-F	26	1.14	9.94	8.72	86.1	
L3153-27	SP2-F-EPH-2	27	1.14	9.94	8.72	86.1	
L3153-28	SP2-G	28	1.13	9.92	8.71	86.2	
L3153-30	SP2-G-EPH-2	29	1.13	9.92	8.71	86.2	
L3153-31	SP2-H	30	1.13	9.94	8.74	86.4	
L3153-33	SP2-H-EPH-2	31	1.13	9.94	8.74	86.4	
L3153-34	SP2-I	32	1.18	9.78	8.7	87.4	
L3153-36	SP2-I-EPH-2	33	1.18	9.78	8.7	87.4	
L3153-37	SP2-J	34	1.12	9.53	8.3	85.4	
L3153-39	SP2-J-EPH-2	35	1.12	9.53	8.3	85.4	
L3153-40	SP2-K	36	1.12	9.9	8.57	84.9	
L3153-42	SP2-K-EPH-2	37	1.12	9.9	8.57	84.9	
L3153-43	SP2-L	38	1.11	9.97	8.71	85.8	
L3153-45	SP2-L-EPH-2	39	1.11	9.97	8.71	85.8	
L3157-01	MH-D19	40	1.13	9.71	7.34	72.4	
L3157-02	MH-D19	41	1.13	9.84	7.18	69.5	
L3157-03	CB-D35	42	1.19	9.51	7.96	81.4	
L3157-04	CB-D35	43	1.12	9.96	8.38	82.1	
L3157-05	CB-D40	44	1.15	9.95	8.69	85.7	
L3157-06	CB-D40	45	1.1	9.72	8.82	89.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$