



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:05
In Date: 07/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 07/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109927

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3184-01	EO-02-070920	37	1.13	9.62	8.95	92.1	
L3184-02	EO-02-070920	38	1.16	9.66	9.06	92.9	
L3184-03	EO-02-070920-1	39	1.14	9.65	9.21	94.8	
L3184-04	EO-02-070920-2	40	1.11	9.86	9.26	93.1	
L3184-05	EO-02-070920-3	41	1.18	9.77	9.07	91.9	
L3250-01	DRUM-DEBRIS-COMP-1	1	1.13	9.56	9.04	93.8	
L3250-03	31087	2	1.14	9.65	9.36	96.6	
L3250-04	31056	3	1.16	9.53	9.27	96.9	
L3250-05	31088	4	1.18	9.83	8.8	88.1	
L3250-06	31090	5	1.16	9.7	8.72	88.5	
L3250-07	31091	6	1.17	9.68	8.59	87.2	
L3250-08	31092	7	1.17	9.67	8.71	88.7	
L3250-09	31133	8	1.13	9.59	8.77	90.3	
L3251-01	BUR-DRUM-COMP	9	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3251-04	31085	10	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3251-05	31131	11	1.17	5.9	5.72	96.2	
L3255-01	CHRT-57	13	1.15	9.96	9.56	95.5	
L3255-02	CHRT-57	14	1.14	9.52	9.08	94.7	
L3255-03	CHRT-33	15	1.14	9.97	3.26	24.0	
L3255-04	CHRT-33	16	1.14	9.97	3.26	24.0	
L3255-05	CHRT-58	17	1.12	9.9	9.34	93.6	
L3255-06	CHRT-58	18	1.18	9.98	9.35	92.8	
L3255-07	RBR-07	19	1.13	9.58	7.78	78.7	
L3255-08	RBR-07	20	1.16	9.86	8.04	79.1	
L3255-09	CHRT-17	21	1.15	9.63	9.26	95.6	
L3255-10	CHRT-17	22	1.18	9.62	9.36	96.9	
L3255-11	CHRT-31	23	1.17	9.78	8.88	89.5	
L3255-12	CHRT-31	24	1.1	9.83	9.22	93.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/10/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:05
In Date: 07/09/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 07/10/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB109927

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3256-01	351	25	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3256-02	337	26	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3256-04	24465	27	1.13	9.67	8.97	91.8	
L3257-01	6FT-MH1	33	1.13	9.66	8.87	90.7	
L3257-02	MH1	34	1.18	9.89	8.65	85.8	
L3257-03	6FT-MH2	35	1.18	5.77	4.13	64.3	
L3257-04	MH2	36	1.19	9.89	7.3	70.2	
L3258-01	20200414-1	12	1.1	9.92	5.63	51.4	
L3259-01	T-2-DRUM	42	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3259-04	C-1-ROLLOFF	43	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3263-01	S5-01-B	28	1.18	9.66	8.52	86.6	
L3263-02	DUP-5	29	1.13	9.65	8.58	87.4	
L3263-03	S5-01-S	30	1.13	9.92	8.91	88.5	
L3263-04	S5-05-B	31	1.12	9.55	8.76	90.6	
L3263-05	S5-05-S	32	1.18	9.66	8.76	89.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$