



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 07/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110025

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3183-06	OK-01-071420	1	1.13	9.73	9.06	92.2	
L3183-07	OK-01-071420	2	1.13	9.82	9.51	96.4	
L3183-08	OK-01-071420-1	3	1.14	9.57	8.62	88.7	
L3183-11	OK-03-071520	80	1.15	9.96	9.51	94.9	
L3183-12	OK-03-071520	81	1.16	9.91	9.22	92.1	
L3183-13	OK-03-071520-1	82	1.19	9.83	9.07	91.2	
L3183-14	OK-03-071520-2	83	1.17	9.97	9.56	95.3	
L3183-15	OK-03-071520-3	84	1.15	9.94	9.16	91.1	
L3183-16	OK-03-071520-4	85	1.14	9.92	9.66	97.0	
L3184-07	EO-01-071520	86	1.18	9.54	8.52	87.8	
L3184-08	EO-01-071520	87	1.13	9.87	8.75	87.2	
L3184-09	EO-01-071520-1	88	1.16	9.76	8.72	87.9	
L3184-10	EO-01-071520-2	89	1.1	9.89	8.89	88.6	
L3184-11	EO-01-071520-3	90	1.19	9.71	8.56	86.5	
L3259-07	C-1-ROLLOFF	38	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3288-01	BS-DRUM	4	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L3289-01	IRON-1	5	1.16	9.93	9.02	89.6	
L3289-02	IRON-1	6	1.13	9.76	8.18	81.7	
L3289-04	IRON-1A	7	1.13	9.76	8.18	81.7	
L3289-05	IRON-1-EPH-2	8	1.13	9.76	8.18	81.7	
L3289-06	IRON-2	9	1.12	9.64	8.22	83.3	
L3289-07	IRON-2	10	1.15	9.96	8.51	83.5	
L3289-09	IRON-2A	11	1.15	9.96	8.51	83.5	
L3289-10	IRON-2-EPH-2	12	1.15	9.96	8.51	83.5	
L3300-01	TP-8A-S1	13	1.12	9.67	8.72	88.9	
L3300-02	TP-8A-S2	14	1.13	9.59	8.56	87.8	
L3300-03	TP-8A-S	15	1.16	9.9	8.92	88.8	
L3300-05	TP-8A-D1	16	1.15	9.94	9.14	90.9	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 07/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110025

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3300-06	TP-8A-D2	17	1.18	9.92	9.11	90.7	
L3300-07	TP-8A-D	18	1.11	9.52	8.68	90.0	
L3301-01	CARLSTADT-SAMPLE-2020	19	1.12	9.88	8.25	81.4	
L3302-01	FK-GF-03-013-A	20	1.16	9.87	8.73	86.9	
L3302-02	FK-GF-03-013-A	21	1.13	9.94	8.78	86.8	
L3302-03	FK-GF-03-013-B	22	1.14	9.95	8.91	88.2	
L3302-04	FK-GF-03-013-B	23	1.16	9.93	8.93	88.6	
L3303-01	FK-GF-03-014-A	24	1.13	9.59	8.48	86.9	
L3303-02	FK-GF-03-014-A	25	1.16	9.78	8.66	87.0	
L3303-03	FK-GF-03-014-B	26	1.18	9.98	8.8	86.6	
L3303-04	FK-GF-03-014-B	27	1.15	9.97	9.12	90.4	
L3304-01	FK-GF-02-058-A	28	1.13	9.94	8.78	86.8	
L3304-02	FK-GF-02-058-A	29	1.13	9.92	8.91	88.5	
L3304-03	FK-GF-02-058-B	30	1.14	9.96	8.9	88.0	
L3304-04	FK-GF-02-058-B	31	1.16	9.53	8.44	87.0	
L3305-01	FK-GF-02-057-A	32	1.17	9.8	8.7	87.3	
L3305-02	FK-GF-02-057-A	33	1.18	9.97	9.06	89.6	
L3305-03	FK-GF-02-057-B	34	1.18	9.96	9.03	89.4	
L3305-04	FK-GF-02-057-B	35	1.18	9.88	8.77	87.2	
L3307-01	ETGI-349	39	1.16	9.74	9.21	93.8	
L3308-01	1	46	1.18	9.88	8.52	84.4	
L3308-02	2	47	1.12	9.58	8.23	84.0	
L3308-03	3	48	1.16	9.56	8.19	83.7	
L3308-04	4	49	1.13	9.53	8.22	84.4	
L3308-05	5	50	1.15	9.88	8.04	78.9	
L3308-06	6	51	1.16	9.6	8.3	84.6	
L3308-07	7	52	1.18	9.77	8.15	81.1	
L3308-08	8	53	1.16	9.9	8.5	84.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 07/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110025

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3308-09	9	54	1.13	9.86	8.35	82.7	
L3308-10	10	55	1.18	9.68	8.15	82.0	
L3308-11	11	56	1.17	9.92	8.56	84.5	
L3308-12	12	57	1.13	9.61	8.04	81.5	
L3308-13	13	58	1.13	9.74	8.29	83.2	
L3308-14	14	59	1.17	9.93	8.37	82.2	
L3308-15	15	60	1.12	9.92	8.6	85.0	
L3308-16	16	61	1.19	9.87	8.25	81.3	
L3308-17	17	62	1.13	9.84	8.47	84.3	
L3308-18	18	63	1.19	9.88	8.25	81.2	
L3308-19	19	64	1.12	9.55	8.24	84.5	
L3308-20	20	65	1.16	9.73	8.13	81.3	
L3308-21	21	66	1.17	9.97	8.64	84.9	
L3308-22	22	67	1.14	9.98	8.56	83.9	
L3308-23	23	68	1.1	9.66	8.14	82.2	
L3308-24	24	69	1.13	9.8	8.36	83.4	
L3310-01	S-1	36	1.13	9.54	8.89	92.3	
L3310-02	S-1A	37	1.16	9.88	9.14	91.5	
L3315-01	TP-16-S1	40	1.15	9.66	8.88	90.8	
L3315-02	TP-16-S2	41	1.16	9.96	9.2	91.4	
L3315-03	TP-16-S	42	1.18	9.76	8.95	90.6	
L3315-05	TP-16-D1	43	1.14	9.94	9.11	90.6	
L3315-06	TP-16-D2	44	1.18	9.97	9.04	89.4	
L3315-07	TP-16-D	45	1.11	9.79	8.85	89.2	
L3316-01	S4-02-B	70	1.1	9.74	8.52	85.9	
L3316-02	S4-03-S	71	1.15	9.95	8.77	86.6	
L3316-03	S4-02-S	72	1.18	9.88	8.71	86.6	
L3316-04	S5-03-B	73	1.13	9.53	8.73	90.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/16/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 07/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110025

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3316-05	S5-03-S	74	1.16	9.98	9.15	90.6	
L3316-06	S2-04-S-01	75	1.13	9.85	8.72	87.0	
L3316-07	S2-06-B-01	76	1.12	9.97	8.07	78.5	
L3316-08	S2-06-S-01	77	1.13	9.7	9.03	92.2	
L3316-09	S2-07-B-01	78	1.17	9.82	8.86	88.9	
L3316-10	S2-08-B-01	79	1.11	9.82	8.89	89.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$