



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 07/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110060

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3188-13	TR-04-071620	1	1.13	9.88	8.99	89.8	
L3188-14	TR-04-071620	2	1.18	9.6	9.06	93.6	
L3188-15	TR-04-071620-1	3	1.17	9.74	9.15	93.1	
L3306-01	WOOSPK-020-0004	4	1.15	9.62	9.03	93.0	
L3306-02	WOOSPK-021-0004	5	1.18	9.94	9.55	95.5	
L3306-03	WOOSPK-025-0004	6	1.13	9.73	8.25	82.8	
L3306-04	WOOSPK-026-0004	7	1.12	9.81	8.69	87.1	
L3306-05	WOOSPK-030-0004	8	1.15	9.76	9.05	91.8	
L3306-06	WOOSPK-031-0004	9	1.16	9.96	9.29	92.4	
L3322-01	31149	10	1.13	9.14	9.02	98.5	
L3322-02	31150	11	1.16	9.79	8.57	85.9	
L3322-03	31151	12	1.14	9.52	9.26	96.9	
L3322-04	31152	13	1.15	9.5	7.53	76.4	
L3322-05	31153	14	1.00	2.00	2.00	100.0	OILY STONE-DEBERIS ,100% SOLIDS
L3322-06	49-53-COMP	15	1.15	9.62	8.22	83.5	
L3323-01	225	16	1.17	9.55	8.35	85.7	
L3323-02	221	17	1.18	9.71	8.49	85.7	
L3324-01	10:1-COMP	18	1.15	9.98	8.56	83.9	
L3324-02	RB47198RT	19	1.16	9.92	8.56	84.5	
L3324-03	CHRT-26948	20	1.16	9.77	8.73	87.9	
L3324-04	CHRT-26948	21	1.18	9.64	8.68	88.7	
L3332-01	S3-04-B-01	22	1.16	9.66	8.46	85.9	
L3332-02	S3-05-B-01	23	1.11	9.75	8.6	86.7	
L3332-03	S4-01-B-01	24	1.19	9.82	8.54	85.2	
L3332-04	S4-01-S-01	25	1.16	9.92	8.62	85.2	
L3335-01	TP-12-S1	26	1.16	9.78	9.1	92.1	
L3335-02	TP-12-S2	27	1.1	9.83	9.05	91.1	
L3335-03	TP-12-S	28	1.14	9.85	9.12	91.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 07/16/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 07/17/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110060

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3335-05	TP-11-S1	29	1.14	9.88	9.1	91.1	
L3335-06	TP-11-S2	30	1.15	9.88	9.1	91.1	
L3335-07	TP-11-S	31	1.13	9.9	9.13	91.2	
L3335-09	TP-11-D1	32	1.16	9.93	8.95	88.8	
L3335-10	TP-11-D2	33	1.19	9.98	8.88	87.5	
L3335-11	TP-11-D	34	1.15	9.75	8.79	88.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$