



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/23/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 14:50
In Date: 07/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 07:50
Out Date: 07/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3349-01	METYC5	1	1.14	9.96	3.46	26.3	
L3349-02	METYC6	2	1.16	9.93	3.00	21.0	
L3349-03	METYC7	3	1.19	9.81	3.04	21.5	
L3349-04	METYC8	4	1.19	9.76	8.18	81.6	
L3349-05	METYC9	5	1.12	9.92	4.94	43.4	
L3349-06	METYD0	6	1.11	9.58	4.63	41.6	
L3349-07	METYD5	7	1.18	9.95	2.63	16.5	
L3349-08	METYD6	8	1.15	9.96	4.29	35.6	
L3349-09	METYD6D	9	1.15	9.96	4.29	35.6	
L3349-10	METYD6S	10	1.15	9.96	4.29	35.6	
L3349-11	METYD7	11	1.13	9.88	5.08	45.1	
L3349-12	METYD8	12	1.14	9.93	5.73	52.2	
L3349-13	METYE1	13	1.16	9.98	2.65	16.9	
L3349-14	METYE2	14	1.14	9.63	3.35	26.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$