



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/23/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 07/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110142

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3369-01	123675	1	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L3369-02	124081	2	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L3369-03	481848	3	1.00	2.00	2.00	100.0	OIL SAMPLE,100% SOLIDS
L3377-01	CON-COMP	27	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3377-02	BRI-COMP	28	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L3377-03	PAINT-COMP	29	1.00	2.00	2.00	100.0	PAINT-CHIPS,SAMPLE,100% SOLIDS
L3378-01	7-22-A	30	1.15	9.98	8.82	86.9	
L3378-02	7-22-B	31	1.16	9.96	7.92	76.8	
L3378-03	7-22-C	32	1.13	9.91	5.67	51.7	
L3378-04	7-22-D	33	1.13	9.93	8.25	80.9	
L3378-05	7-22-E	34	1.18	9.81	7.75	76.1	
L3378-06	7-22-F	35	1.15	9.96	8.00	77.8	
L3378-07	7-22-G	36	1.17	9.7	7.64	75.8	
L3378-08	7-22-H	37	1.18	9.97	8.05	78.2	
L3378-09	7-22-I	38	1.17	9.98	5.68	51.2	
L3378-10	7-22-J	39	1.14	9.98	8.81	86.8	
L3378-11	7-22-K	40	1.13	9.9	7.00	66.9	
L3378-12	7-22-L	41	1.14	9.96	8.00	77.8	
L3378-13	DRUM-COMP	42	1.17	9.94	8.23	80.5	
L3378-14	DRUM-COMP	43	1.17	9.93	7.44	71.6	
L3379-01	ETGI-361	4	1.14	9.98	8.45	82.7	
L3380-01	VNJ-221	5	1.17	9.95	9.33	92.9	
L3380-03	VNJ-241	6	1.17	9.84	8.18	80.9	
L3380-05	VNJ-242	7	1.13	9.86	8.75	87.3	
L3381-01	357-GRID-A4-A6	8	1.19	9.72	7.6	75.1	
L3381-02	357-GRID-A4-A6	9	1.18	9.91	8.79	87.2	
L3381-04	347-GRID-B1-B2	10	1.15	9.52	8.65	89.6	
L3381-05	347-GRID-B1-B2	11	1.14	9.94	9.16	91.1	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/23/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 07/22/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/23/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110142

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3381-07	364-GRID-A7-A8	12	1.17	9.98	8.64	84.8	
L3381-08	364-GRID-A7-A8	13	1.19	9.75	8.03	79.9	
L3381-10	359-GRID-A1-A3	14	1.19	9.96	9.02	89.3	
L3381-11	359-GRID-A1-A3	15	1.1	9.77	8.85	89.4	
L3381-13	339-GRID-B3-B4	16	1.13	9.68	9.17	94.0	
L3381-14	339-GRID-B3-B4	17	1.19	9.94	9.45	94.4	
L3381-16	365-GRID-B6-B8	18	1.15	9.6	8.75	89.9	
L3381-17	365-GRID-B6-B8	19	1.18	9.82	9.04	91.0	
L3381-19	VNJ-225-GRID-B5	20	1.18	9.92	9.28	92.7	
L3381-20	VNJ-225-GRID-B5	21	1.1	9.81	9.14	92.3	
L3403-01	T-1	22	1.11	9.98	8.97	88.6	
L3403-02	T-2	23	1.17	9.9	8.8	87.4	
L3403-03	T-3	24	1.13	9.93	8.86	87.8	
L3403-04	T-4	25	1.15	9.97	8.95	88.4	
L3403-05	T-5	26	1.17	9.97	8.65	85.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$