



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 07/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110184

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3351-06	SAMPLE-1A	6	1.18	9.78	9.02	91.2	
L3351-07	SAMPLE-1B	7	1.19	9.63	9.23	95.3	
L3428-01	TP-1-S1	1	1.14	9.83	9.08	91.4	
L3428-02	TP-1-S2	2	1.13	9.8	8.98	90.5	
L3428-03	TP-1-S	3	1.14	9.96	9.2	91.4	
L3430-01	TS1	4	1.18	9.66	9.08	93.2	
L3430-02	TS2	5	1.12	9.52	8.45	87.3	
L3432-01	31164	71	1.12	9.94	9.52	95.2	
L3432-02	31165-S	72	1.15	9.96	8.13	79.2	
L3432-04	31166	73	1.18	9.56	9.25	96.3	
L3432-05	31167	74	1.17	9.92	9.48	95.0	
L3432-06	64-66-67	75	1.13	9.96	9.59	95.8	
L3434-01	11986	76	1.18	9.94	9.67	96.9	
L3434-03	11936	77	1.15	9.98	6.55	61.2	
L3435-01	SP-1	78	1.13	9.9	8.82	87.7	
L3435-02	SP-1	79	1.17	9.8	8.78	88.2	
L3435-03	SP-2	80	1.13	9.9	8.74	86.8	
L3435-04	SP-2	81	1.15	9.96	9.46	94.3	
L3435-05	SP-3	82	1.1	9.74	8.85	89.7	
L3435-06	SP-3	83	1.19	9.89	9.00	89.8	
L3436-01	1	15	1.18	9.76	9.09	92.2	
L3436-02	2	16	1.17	9.77	8.97	90.7	
L3436-03	3	17	1.13	9.54	8.73	90.4	
L3436-04	4	18	1.18	9.79	9.06	91.5	
L3436-05	5	19	1.15	9.86	9.00	90.1	
L3436-06	6	20	1.15	9.75	8.97	90.9	
L3436-07	7	21	1.19	9.93	9.22	91.9	
L3436-08	8	22	1.17	9.97	9.3	92.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 07/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110184

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3436-09	9	23	1.18	9.83	9.15	92.1	
L3436-10	10	24	1.14	9.92	9.16	91.3	
L3436-11	11	25	1.18	9.78	9.06	91.6	
L3436-12	12	26	1.18	9.93	9.27	92.5	
L3436-13	13	27	1.17	9.6	8.82	90.7	
L3436-14	14	28	1.12	9.94	9.15	91.0	
L3436-15	15	29	1.14	9.91	9.2	91.9	
L3436-16	16	30	1.16	9.59	8.89	91.7	
L3436-17	17	31	1.18	9.66	8.82	90.1	
L3436-18	18	32	1.17	9.86	8.72	86.9	
L3436-19	19	33	1.17	9.86	9.18	92.2	
L3436-20	20	34	1.16	9.98	9.16	90.7	
L3436-21	21	35	1.18	9.77	9.05	91.6	
L3436-22	22	36	1.15	9.95	9.18	91.3	
L3439-01	20200438-3001A	8	1.15	9.98	5.75	52.1	
L3439-02	20200439-3001B	9	1.13	9.71	5.52	51.2	
L3439-03	20200440-3001C	10	1.18	9.6	5.51	51.4	
L3439-04	20200441-3002A	11	1.19	9.68	5.54	51.2	
L3439-05	20200442-3002B	12	1.15	9.94	5.72	52.0	
L3439-06	20200443-3002C	13	1.15	9.96	5.72	51.9	
L3439-07	20200444-3002D	14	1.13	9.55	5.46	51.4	
L3440-05	SVOC-GPC-BLANK	37	1.00	2.00	2.00	100.0	
L3440-06	PEST-GPC-BLANK	38	1.00	2.00	2.00	100.0	
L3440-07	PEST-GPC-BLANK-SPIKE	39	1.00	2.00	2.00	100.0	
L3440-08	PCB-GPC-BLANK	40	1.00	2.00	2.00	100.0	
L3440-09	PCB-GPC-BLANK-SPIKE	41	1.00	2.00	2.00	100.0	
L3440-10	SVOC-GPC2-BLANK	42	1.00	2.00	2.00	100.0	
L3440-11	PEST-GPC2-BLANK	43	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 07/24/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/25/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110184

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3440-12	PEST-GPC2-BLANK-SPIKE	44	1.00	2.00	2.00	100.0	
L3440-13	PCB-GPC2-BLANK	45	1.00	2.00	2.00	100.0	
L3440-14	PCB-GPC2-BLANK-SPIKE	46	1.00	2.00	2.00	100.0	
L3444-01	P001-CC44-B1-1824-01	57	1.11	9.63	9.17	94.6	
L3444-02	P001-CC46-B1-1218-01	58	1.15	9.97	8.17	79.6	
L3444-03	P001-CH45-W1-1824-02	59	1.14	9.96	9.32	92.7	
L3444-04	P001-CH45-W1-1824-01	60	1.13	9.66	9.06	93.0	
L3444-05	P001-CH45-W1-1824-01MS	61	1.13	9.66	9.06	93.0	
L3444-06	P001-CH45-W1-1824-01MSD	62	1.13	9.66	9.06	93.0	
L3444-07	P001-CI45-W1-1824-01	63	1.17	9.94	8.04	78.3	
L3445-01	P001-CC44-B1-1824-01	64	1.11	9.63	9.17	94.6	
L3445-02	P001-CC46-B1-1218-01	65	1.15	9.97	8.17	79.6	
L3445-03	P001-CH45-W1-1824-02	66	1.14	9.96	9.32	92.7	
L3445-04	P001-CH45-W1-1824-01	67	1.13	9.66	9.06	93.0	
L3445-05	P001-CH45-W1-1824-01D	68	1.13	9.66	9.06	93.0	
L3445-06	P001-CH45-W1-1824-01S	69	1.13	9.66	9.06	93.0	
L3445-07	P001-CI45-W1-1824-01	70	1.17	9.94	8.04	78.3	
L3449-01	FK-GF-02-059-A	47	1.14	9.81	8.92	89.7	
L3449-02	FK-GF-02-059-A	48	1.15	9.96	9.04	89.6	
L3449-03	FK-GF-02-059-B	49	1.17	9.75	8.72	88.0	
L3449-04	FK-GF-02-059-B	50	1.19	9.81	8.88	89.2	
L3451-01	FK-GF-03-015-A	51	1.14	9.97	9.09	90.0	
L3451-02	FK-GF-03-015-A	52	1.14	9.98	9.03	89.3	
L3451-03	FK-GF-03-015-B	53	1.13	9.68	8.83	90.1	
L3451-04	FK-GF-03-015-B	54	1.16	9.82	8.91	89.5	
L3453-01	FK-PS-01-067	55	1.11	9.76	8.4	84.3	
L3453-02	FK-PS-01-067	56	1.17	9.86	8.43	83.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$