



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/28/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 14:50
In Date: 07/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 07/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110215

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3458-01	GAY29	1	1.15	9.95	8.52	83.8	
L3458-02	GAY30	2	1.13	9.71	8.09	81.1	
L3458-03	GAY31	3	1.16	9.97	8.09	78.7	
L3458-04	GAY32	4	1.14	9.94	8.82	87.3	
L3458-05	GAY33	5	1.17	9.98	7.69	74.0	
L3458-06	GAY34	6	1.15	9.9	7.77	75.7	
L3458-07	GAY35	7	1.12	9.92	8.83	87.6	
L3458-08	GAY36	8	1.14	9.7	7.86	78.5	
L3458-09	GAY36MS	9	1.14	9.7	7.86	78.5	
L3458-10	GAY36MSD	10	1.14	9.7	7.86	78.5	
L3458-11	GAY37	11	1.18	9.66	7.88	79.0	
L3458-12	GAY38	12	1.16	9.82	8.64	86.4	
L3458-13	GAY39	13	1.14	9.95	8.5	83.5	
L3458-14	GAY40	14	1.16	9.85	8.65	86.2	
L3458-15	GAY41	15	1.11	9.65	8.67	88.5	
L3458-16	GAY42	16	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3458-17	VHBLK01	17	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$