



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 07/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:02
Out Date: 07/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110270

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3034-22	TP-3	1	1.18	9.91	8.85	87.9	
L3180-01	HD-02-072920	2	1.14	9.9	8.77	87.1	
L3180-02	HD-02-072920	3	1.18	9.92	8.88	88.1	
L3180-03	HD-02-072920-1	4	1.15	9.77	8.73	87.9	
L3478-01	ROLLOFF-VNJ-213	35	1.15	9.57	8.7	89.7	
L3478-02	ROLLOFF-VNJ-213	36	1.18	9.55	7.56	76.2	
L3479-01	REEL-KMH272B-END-1	37	1.00	2.00	2.00	100.0	
L3479-02	REEL-KMH272B-END-2	38	1.00	2.00	2.00	100.0	
L3479-03	REEL-HEF239Y-END-1	39	1.00	2.00	2.00	100.0	
L3479-04	REEL-HEF239Y-END-2	40	1.00	2.00	2.00	100.0	
L3480-01	X068-B8	5	1.19	9.75	9.03	91.6	
L3480-02	X068-B9	6	1.15	9.56	8.02	81.7	
L3480-03	X068-B5	7	1.13	9.75	8.00	79.7	
L3481-01	ROLLOFF-72-12012	41	1.17	9.61	8.78	90.2	
L3481-02	ROLLOFF-72-12012	42	1.15	9.95	8.69	85.7	
L3483-01	S-1	30	1.18	9.7	9.13	93.3	
L3483-03	20C023	31	1.00	2.00	2.00	100.0	concrete SAMPLE,100% SOLIDS
L3484-01	30982	32	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3484-02	31156	33	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3484-03	82-56	34	1.00	2.00	2.00	100.0	OILY DEBERIS,100% SOLIDS
L3489-01	CS-1	8	1.13	9.53	9.34	97.7	
L3489-02	CS-2	9	1.14	9.62	9.44	97.9	
L3489-03	CS-3	10	1.18	9.64	9.43	97.5	
L3489-04	CS-4	11	1.17	9.79	9.61	97.9	
L3489-05	CS-5	12	1.16	9.56	9.33	97.3	
L3489-06	CS-6	13	1.18	9.59	9.4	97.7	
L3489-07	CS-7	14	1.14	9.63	9.45	97.9	
L3489-08	CS-8	15	1.15	9.75	9.52	97.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 07/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:02
Out Date: 07/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110270

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3489-09	CS-9	16	1.18	9.57	9.4	98.0	
L3489-10	CS-10	17	1.19	9.63	9.48	98.2	
L3489-11	CS-11	18	1.15	9.69	9.48	97.5	
L3489-12	CS-12	19	1.12	9.72	9.52	97.7	
L3491-01	TP-1-S1	20	1.19	9.6	8.47	86.6	
L3491-02	TP-1-S2	21	1.14	9.95	8.76	86.5	
L3491-03	TP-1-S	22	1.17	9.54	8.58	88.5	
L3491-05	TP-1-D1	23	1.16	9.88	8.46	83.7	
L3491-06	TP-1-D2	24	1.19	9.8	8.53	85.2	
L3491-07	TP-1-D3	25	1.16	9.97	9.47	94.3	
L3491-08	TP-1-D	26	1.12	9.74	8.55	86.2	
L3491-10	TP-3-S1	27	1.14	9.62	8.54	87.3	
L3491-11	TP-3-S2	28	1.11	9.58	8.17	83.4	
L3491-12	TP-3-S	29	1.15	9.55	8.64	89.2	
L3494-01	JTY-SB-05-2-3	43	1.18	9.69	9.08	92.8	
L3494-02	JTY-SB-105-2-3	44	1.13	9.82	9.03	90.9	
L3494-03	JTY-SB-03-2-3	45	1.12	9.54	8.79	91.1	
L3494-04	JTY-SB-02-2-3	46	1.17	9.96	9.18	91.1	
L3494-05	JTY-SB-04-2-3	47	1.13	9.89	9.11	91.1	
L3494-06	JTY-SB-01-2-3	48	1.14	9.51	8.8	91.5	
L3494-07	JTY-SB-02-9-10	49	1.17	9.94	9.15	91.0	
L3494-08	JTY-SB-03-9-10	50	1.18	8.25	7.64	91.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$