



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/31/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 07/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 07/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110292

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3482-01	EFB082M-END-1	15	1.00	2.00	2.00	100.0	
L3482-02	EFB082M-END-2	16	1.00	2.00	2.00	100.0	
L3482-03	KMA903F-END-1	17	1.00	2.00	2.00	100.0	
L3482-04	KMA903F-END-2	18	1.00	2.00	2.00	100.0	
L3482-05	KMH408K-END-1	19	1.00	2.00	2.00	100.0	
L3482-06	KMH408K-END-2	20	1.00	2.00	2.00	100.0	
L3482-07	JEB403X-END-1	21	1.00	2.00	2.00	100.0	
L3482-08	JEB403X-END-2	22	1.00	2.00	2.00	100.0	
L3503-01	20200449-2+90	1	1.18	9.9	8.44	83.3	
L3507-01	SU-01-073020	9	1.18	9.66	8.77	89.5	
L3507-02	SU-01-073020	10	1.15	9.92	8.82	87.5	
L3507-03	SU-01-073020-1	11	1.14	9.96	9.36	93.2	
L3507-04	SU-01-073020-2	12	1.13	9.97	9.44	94.0	
L3507-05	SU-01-073020-3	13	1.15	9.5	8.73	90.8	
L3507-06	SU-01-073020-4	14	1.13	9.75	9.52	97.3	
L3511-02	TP-1	23	1.18	9.68	9.14	93.6	
L3511-05	TP-1-EPH-2	24	1.18	9.82	9.02	90.7	
L3511-07	TP-2	25	1.18	9.56	8.92	92.4	
L3511-10	TP-2-EPH-2	26	1.1	9.86	9.06	90.9	
L3511-12	TP-3	27	1.14	9.88	9.14	91.5	
L3511-15	TP-3-EPH-2	28	1.15	9.94	9.26	92.3	
L3511-17	SP	29	1.14	9.96	9.47	94.4	
L3511-20	SP-EPH-2	30	1.14	9.83	9.54	96.7	
L3513-01	X068-B2	2	1.15	9.75	9.06	92.0	
L3513-02	X068-S2	3	1.19	9.67	9.22	94.7	
L3513-03	X068-DUP2	4	1.19	9.66	8.06	81.1	
L3513-04	X068-B1	5	1.15	9.95	8.95	88.6	
L3513-05	X068-S1	6	1.17	9.96	9.04	89.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/31/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 07/30/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 07/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110292

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3513-06	X068-S5	7	1.13	9.57	8.09	82.5	
L3513-07	X068-DUP1	8	1.16	9.63	8.93	91.7	
L3518-02	TP-4	31	1.19	9.74	9.14	93.0	
L3518-05	TP-4-EPH-2	32	1.15	9.96	9.61	96.0	
L3518-07	TP-5	33	1.15	9.96	9.57	95.6	
L3518-10	TP-5-EPH-2	34	1.16	9.93	9.47	94.8	
L3518-12	TP-6	35	1.16	9.79	9.16	92.7	
L3518-15	TP-6-EPH-2	36	1.1	9.9	8.42	83.2	
L3518-17	TP-7	37	1.16	9.96	8.69	85.6	
L3518-20	TP-7-EPH-2	38	1.15	9.98	8.69	85.4	
L3518-22	TP-8	39	1.15	9.94	9.55	95.6	
L3518-25	TP-8-EPH-2	40	1.1	9.92	9.2	91.8	
L3518-27	TP-9	41	1.13	9.9	9.61	96.7	
L3518-30	TP-9-EPH-2	42	1.13	9.91	9.33	93.4	
L3518-32	TP-10	43	1.15	9.94	8.98	89.1	
L3518-35	TP-10-EPH-2	44	1.17	9.96	9.22	91.6	
L3518-37	TP-11	45	1.12	9.75	9.27	94.4	
L3518-40	TP-11-EPH-2	46	1.14	9.97	9.55	95.2	
L3518-42	TP-12	47	1.18	9.89	8.94	89.1	
L3518-45	TP-12-EPH-2	48	1.18	9.73	8.73	88.3	
L3518-47	TP-13	49	1.14	9.94	8.78	86.8	
L3518-50	TP-13-EPH-2	50	1.13	9.62	8.86	91.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$