



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 8/11/2020

OVENTEMP IN Celsius(°C): 108

Time IN: 16:40

In Date: 08/10/2020

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:00

Out Date: 08/11/2020

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB110461

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3609-04	X068-B15-02	1	1.17	9.97	8.06	78.3	
L3623-02	SB-2	2	1.15	9.81	8.05	79.7	
L3623-03	STP-7A	3	1.19	9.76	9.26	94.2	
L3624-02	SS-TP-1	8	1.18	9.74	9.1	92.5	
L3624-05	SS-TP-1-EPH-2	9	1.13	9.75	9.06	92.0	
L3624-07	SS-TP-2	10	1.15	9.96	9.7	97.0	
L3624-10	SS-TP-2-EPH-2	11	1.12	9.92	9.39	94.0	
L3624-12	SS-TP-3	12	1.18	9.66	9.17	94.2	
L3624-15	SS-TP-3-EPH-2	13	1.14	9.96	9.27	92.2	
L3624-17	SS-TP-4	14	1.1	9.83	8.82	88.4	
L3624-20	SS-TP-4-EPH-2	15	1.13	9.72	8.78	89.1	
L3625-01	PN-SP	16	1.16	9.79	9.08	91.8	
L3625-02	PN-SP-EPH-2	17	1.17	9.8	9.13	92.2	
L3627-01	WH-CON-A	18	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L3627-02	WH-CON-B	19	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
L3638-01	X445-S1	4	1.12	9.58	8.34	85.3	
L3638-02	X445-S2	5	1.14	9.84	8.79	87.9	
L3638-03	X445-B1	6	1.18	9.66	9.15	94.0	
L3638-04	X445-DUP-1	7	1.16	9.67	8.67	88.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$