



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/12/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 08/11/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 08/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110467

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3632-01	BFM79	1	1.13	9.58	8.63	88.8	
L3632-02	BFM80	2	1.14	9.78	9.32	94.7	
L3632-03	BFM81	3	1.15	9.9	9.48	95.2	
L3632-04	BFM82	4	1.15	9.98	9.51	94.7	
L3632-05	BFM83	5	1.14	9.88	9.47	95.3	
L3632-06	BFM84	6	1.12	9.86	9.34	94.1	
L3632-07	BFM85	7	1.18	9.92	8.51	83.9	
L3632-08	BFM86	8	1.19	9.83	9.16	92.2	
L3632-09	BFM87	9	1.1	9.63	8.96	92.1	
L3632-10	BFM88	10	1.13	9.69	8.69	88.3	
L3632-11	BFM89	11	1.15	9.88	8.9	88.8	
L3632-12	BFM90	12	1.15	9.87	8.7	86.6	
L3632-13	BFM91	13	1.16	9.62	9.02	92.9	
L3632-14	BFM92	14	1.17	9.8	9.26	93.7	
L3632-15	BFM93	15	1.18	9.64	9.16	94.3	
L3632-16	BFM94	16	1.13	9.79	9.22	93.4	
L3632-17	BFM95	17	1.18	9.97	8.77	86.3	
L3632-18	BFM96	18	1.14	9.5	8.75	91.0	
L3632-19	BFM97	19	1.16	9.96	9.47	94.4	
L3632-20	BFM98	20	1.17	9.98	9.52	94.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$