



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/13/2020

OVENTEMP IN Celsius(°C): 108
Time IN: 14:35
In Date: 08/12/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 08/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110490

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3633-01	BFM99	1	1.14	9.88	9.28	93.1	
L3633-02	BFMA0	2	1.15	9.66	8.55	87.0	
L3633-03	BFMA1	3	1.17	9.77	9.15	92.8	
L3633-04	BFMA2	4	1.12	9.53	8.97	93.3	
L3633-05	BFMA3	5	1.15	9.68	9.13	93.6	
L3633-06	BFMA4	6	1.17	9.94	9.36	93.4	
L3633-07	BFMA5	7	1.14	9.94	9.55	95.6	
L3633-08	BFMA6	8	1.14	9.65	8.8	90.0	
L3633-09	BFMA7	9	1.18	9.95	9.46	94.4	
L3633-10	BFMD8	10	1.15	9.87	9.45	95.2	
L3633-11	BFMD9	11	1.12	9.57	9.19	95.5	
L3633-12	BFME0	12	1.17	9.71	9.22	94.3	
L3633-13	BFME1	13	1.13	9.66	8.46	85.9	
L3633-14	BFME2	14	1.15	9.75	9.03	91.6	
L3633-15	BFME3	15	1.17	9.96	9.35	93.1	
L3633-16	BFME4	16	1.16	9.71	9.09	92.7	
L3633-17	BFME5	17	1.12	9.94	9.22	91.8	
L3633-18	BFME6	18	1.14	9.55	8.44	86.8	
L3633-19	BFME7	19	1.1	9.82	8.58	85.8	
L3633-20	BFME8	20	1.13	9.61	9.05	93.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$