



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/14/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/13/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110519

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3686-13	MX2N03	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-14	MX2N03D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-15	MX2N03S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-16	MX2N04	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-17	MX2N04D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-18	MX2N04S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-19	MX2N05	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-20	MX2N05D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-21	MX2N05S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-22	MX2N06	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-23	MX2N06D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3686-24	MX2N06S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$