



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/17/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/14/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 08/15/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB110534

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3673-01	MGAYD0	1	1.15	9.59	9.32	96.8	
L3673-02	MGAYD1	2	1.18	9.87	9.42	94.8	
L3673-03	MGAYD2	3	1.17	9.67	9.37	96.5	
L3673-04	MGAYD3	4	1.18	9.66	9.47	97.8	
L3673-05	MGAYD4	5	1.12	9.52	9.33	97.7	
L3673-06	MGAYD5	6	1.14	9.72	9.44	96.7	
L3673-07	MGAYD6	7	1.17	9.63	9.4	97.3	
L3673-08	MGAYD7	8	1.14	9.88	9.54	96.1	
L3673-09	MGAYD8	9	1.15	9.55	9.32	97.3	
L3673-10	MGAYD9	10	1.17	9.9	9.69	97.6	
L3673-11	MGAYD9D	11	1.17	9.9	9.69	97.6	
L3673-12	MGAYD9S	12	1.17	9.9	9.69	97.6	
L3673-13	MGAYE0	13	1.18	9.73	9.4	96.1	
L3673-14	MGAYE1	14	1.17	9.86	9.24	92.9	
L3673-15	MGAYE2	15	1.1	9.53	9.29	97.2	
L3673-16	MGAYE3	16	1.15	9.88	9.63	97.1	
L3673-17	MGAYE4	17	1.17	9.91	9.7	97.6	
L3673-18	MGAYE5	18	1.00	2.00	2.00	100.0	P.E. SAMPLE
L3673-19	MGAYE6	19	1.12	9.69	9.5	97.8	
L3673-20	MGAYE7	20	1.19	9.89	9.55	96.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$