



PERCENT SOLID

Supervisor: apatel  
Analyst: JIGNESH  
Date: 8/18/2020

OVENTEMP IN Celsius(°C): 107  
Time IN: 17:20  
In Date: 08/17/2020  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:15  
Out Date: 08/18/2020  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %SOLIDS-OVEN

QC:LB110563

| Lab ID   | Client SampleID           | Dish # | Dish Wt(g) (A) | Dish + Sample Wt(g) (B) | Dish+Dry Sample Wt(g) (C) | % Solid | Comments    |
|----------|---------------------------|--------|----------------|-------------------------|---------------------------|---------|-------------|
| L3496-01 | AU-06-081720              | 53     | 1.15           | 9.96                    | 9.14                      | 90.7    |             |
| L3496-02 | AU-06-081720              | 54     | 1.13           | 9.53                    | 8.41                      | 86.7    |             |
| L3496-03 | AU-06-081720-1            | 55     | 1.14           | 9.81                    | 8.87                      | 89.2    |             |
| L3496-04 | AU-06-081720-2            | 56     | 1.13           | 9.77                    | 9.00                      | 91.1    |             |
| L3496-05 | AU-06-081720-3            | 57     | 1.14           | 9.96                    | 9.15                      | 90.8    |             |
| L3689-01 | TP-1-4-8FT                | 1      | 1.15           | 9.97                    | 9.15                      | 90.7    |             |
| L3689-02 | TP-1-4-8FT                | 2      | 1.16           | 9.61                    | 8.58                      | 87.8    |             |
| L3689-04 | TP-5                      | 3      | 1.18           | 9.77                    | 8.49                      | 85.1    |             |
| L3689-05 | TP-5                      | 4      | 1.15           | 9.97                    | 9.22                      | 91.5    |             |
| L3689-07 | TP-4                      | 5      | 1.14           | 9.9                     | 9.2                       | 92.0    |             |
| L3689-08 | TP-4                      | 6      | 1.15           | 9.93                    | 9.18                      | 91.5    |             |
| L3689-10 | TP-3                      | 7      | 1.18           | 9.89                    | 8.15                      | 80.0    |             |
| L3689-11 | TP-3                      | 8      | 1.17           | 9.79                    | 8.7                       | 87.4    |             |
| L3700-01 | TB-01-081420              | 9      | 1.00           | 2.00                    | 2.00                      | 100.0   | T.B.        |
| L3700-02 | S406-K-SO-1-1.5-081420    | 10     | 1.15           | 9.98                    | 8.82                      | 86.9    |             |
| L3700-03 | L3700-01MS                | 11     | 1.15           | 9.98                    | 8.82                      | 86.9    |             |
| L3700-04 | L3700-01MSD               | 12     | 1.15           | 9.98                    | 8.82                      | 86.9    |             |
| L3700-05 | S405-K-SO-1-1.5-081420    | 13     | 1.13           | 9.74                    | 8.18                      | 81.9    |             |
| L3700-06 | S407-K-SO-1-1.5-081420    | 14     | 1.14           | 9.68                    | 7.48                      | 74.2    |             |
| L3700-07 | S407-K-SO-1-1.5-081420-FD | 15     | 1.16           | 9.92                    | 7.53                      | 72.7    |             |
| L3700-08 | S408-K-SO-1-1.5-081420    | 16     | 1.16           | 9.75                    | 8.59                      | 86.5    |             |
| L3700-09 | S709-N-SO-0-0.5-081420    | 17     | 1.15           | 9.96                    | 9.03                      | 89.4    |             |
| L3704-01 | TR1-PLYWOOD-TOP           | 18     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |
| L3704-02 | TR1-PLYWOOD-BOT           | 19     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |
| L3704-03 | TR2-PLYWOOD-TOP           | 20     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |
| L3704-04 | TR2-PLYWOOD-BOT           | 21     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |
| L3705-01 | TR1-WIPE-1                | 22     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |
| L3705-02 | TR1-WIPE-2                | 23     | 1.00           | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE |

**PERCENT SOLID**

**Supervisor:** apatel  
**Analyst:** JIGNESH  
**Date:** 8/18/2020

**OVENTEMP IN Celsius(°C):** 107  
**Time IN:** 17:20  
**In Date:** 08/17/2020  
**Weight Check 1.0g:** 1.00  
**Weight Check 10g:** 10.00  
**OvenID:** M OVEN-1

**OVENTEMP OUT Celsius(°C):** 103  
**Time OUT:** 08:15  
**Out Date:** 08/18/2020  
**Weight Check 1.0g:** 1.00  
**Weight Check 10g:** 10.00  
**BalanceID:** M SC-4  
**Thermometer ID:** %SOLIDS-OVEN

QC:LB110563

| Lab ID   | Client SampleID               | Dish # | Dish Wt (g) (A) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments      |
|----------|-------------------------------|--------|-----------------|--------------------------|----------------------------|---------|---------------|
| L3705-03 | TR1-WIPE-3                    | 24     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3705-04 | TR1-WIPE-4                    | 25     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3705-05 | TR2-WIPE-1                    | 26     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3705-06 | TR2-WIPE-2                    | 27     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3705-07 | TR2-WIPE-3                    | 28     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3705-08 | TR2-WIPE-4                    | 29     | 1.00            | 2.00                     | 2.00                       | 100.0   | WIPE SAMPLE   |
| L3711-01 | TT-025-IDWSO-20200814         | 30     | 1.14            | 9.63                     | 8.03                       | 81.2    |               |
| L3711-02 | TT-026-IDWSO-20200814         | 31     | 1.15            | 9.98                     | 8.48                       | 83.0    |               |
| L3712-01 | CHRT-27137                    | 33     | 1.18            | 9.77                     | 7.96                       | 78.9    |               |
| L3712-02 | CHRT-27137                    | 34     | 1.13            | 9.59                     | 7.76                       | 78.4    |               |
| L3712-03 | CHRT-26572                    | 35     | 1.19            | 9.88                     | 9.3                        | 93.3    |               |
| L3712-04 | CHRT-26572                    | 36     | 1.19            | 9.74                     | 9.22                       | 93.9    |               |
| L3712-05 | VB-25203                      | 37     | 1.18            | 9.61                     | 8.26                       | 84.0    |               |
| L3712-06 | VB-25203                      | 38     | 1.16            | 9.84                     | 8.48                       | 84.3    |               |
| L3712-07 | VB-27207                      | 39     | 1.12            | 9.7                      | 8.12                       | 81.6    |               |
| L3712-08 | VB-27207                      | 40     | 1.17            | 9.9                      | 7.68                       | 74.6    |               |
| L3712-09 | V-23615                       | 41     | 1.1             | 9.67                     | 8.36                       | 84.7    |               |
| L3712-10 | V-23615                       | 42     | 1.16            | 9.77                     | 8.46                       | 84.8    |               |
| L3712-11 | CHRT-25433                    | 43     | 1.19            | 9.7                      | 7.31                       | 71.9    |               |
| L3712-12 | CHRT-25433                    | 44     | 1.16            | 9.66                     | 7.94                       | 79.8    |               |
| L3712-13 | CHRT-26982                    | 45     | 1.19            | 9.97                     | 7.17                       | 68.1    |               |
| L3712-14 | CHRT-26982                    | 46     | 1.17            | 9.78                     | 6.72                       | 64.5    |               |
| L3717-01 | BP-VPB179-GW-738-740          | 32     | 1.18            | 9.98                     | 2.86                       | 19.1    | SLUDGE SAMPLE |
| L3719-01 | 20200506-DISCRETE-004         | 47     | 1.15            | 9.95                     | 5.91                       | 54.1    |               |
| L3719-02 | 20200507-DISCRETE-005         | 48     | 1.15            | 9.91                     | 4.32                       | 36.2    |               |
| L3719-03 | 20200509-COMPOSITE-A          | 49     | 1.17            | 9.97                     | 4.04                       | 32.6    |               |
| L3719-04 | 20200510-AMENDED-COMP-A       | 50     | 1.16            | 9.7                      | 4.36                       | 37.5    |               |
| L3719-05 | 20200511-AMENDED-DISCRETE-004 | 51     | 1.12            | 9.83                     | 5.9                        | 54.9    |               |



PERCENT SOLID

Supervisor: apatel  
Analyst: JIGNESH  
Date: 8/18/2020

OVENTEMP IN Celsius(°C): 107  
Time IN: 17:20  
In Date: 08/17/2020  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:15  
Out Date: 08/18/2020  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %SOLIDS-OVEN

QC:LB110563

| Lab ID   | Client SampleID               | Dish # | Dish Wt (g) (A) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments |
|----------|-------------------------------|--------|-----------------|--------------------------|----------------------------|---------|----------|
| L3719-06 | 20200512-AMENDED-DISCRETE-005 | 52     | 1.17            | 9.69                     | 4.71                       | 41.5    |          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$