



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/20/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110605

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3499-01	HD-01-081820	1	1.15	9.91	9.25	92.5	
L3499-02	HD-01-081820	2	1.16	9.83	8.95	89.9	
L3499-03	HD-01-081820-1	3	1.14	9.75	9.37	95.6	
L3502-11	OK-02-081920	34	1.15	9.83	9.04	90.9	
L3502-12	OK-02-081920	35	1.15	9.88	9.22	92.4	
L3502-13	OK-02-081920-1	36	1.13	9.78	9.06	91.7	
L3502-14	OK-02-081920-2	37	1.16	9.53	8.89	92.4	
L3504-07	EO-01-081920	38	1.14	9.55	8.74	90.4	
L3504-08	EO-01-081920	39	1.15	9.94	9.14	90.9	
L3504-09	EO-01-081920-1	40	1.18	9.66	8.9	91.0	
L3504-10	EO-01-081920-2	41	1.1	9.72	8.86	90.0	
L3504-11	EO-01-081920-3	42	1.16	9.97	9.08	89.9	
L3715-01	226	4	1.00	2.00	2.00	100.0	OIL SAMPLE
L3720-03	MP-081920-B10-0-2-COMP	43	1.18	9.54	8.93	92.7	
L3720-05	MP-081920-B11-0-2-COMP	44	1.19	9.94	8.9	88.1	
L3720-07	MP-081920-BD-0-2-COMP	45	1.11	9.93	8.41	82.8	
L3720-09	MP-081920-B9-0-2-COMP	46	1.19	9.85	8.26	81.6	
L3729-01	SB-3- (0-1)	5	1.18	9.63	9.17	94.6	
L3729-02	SB-4- (3-4)	6	1.15	9.96	9.25	91.9	
L3729-03	SB-5- (4-5)	7	1.17	9.92	9.44	94.5	
L3729-04	SB-2- (3-4)	8	1.16	9.6	8.93	92.1	
L3729-05	SB-1- (1-2)	9	1.18	9.77	8.75	88.1	
L3729-06	SB-5- (3-4)	10	1.18	9.98	9.55	95.1	
L3729-07	SB-2- (4-5)	11	1.13	9.85	8.41	83.5	
L3729-08	SB-3- (4-5)	12	1.14	9.92	9.15	91.2	
L3729-09	SB-2- (1-2)	13	1.18	9.71	9.43	96.7	
L3729-10	SB-4- (1-2)	14	1.16	9.58	9.00	93.1	
L3729-11	SB-5- (2-3)	15	1.17	9.52	9.01	93.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/20/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/19/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/20/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110605

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3729-12	SB-3-(2-3)	16	1.12	9.88	9.48	95.4	
L3729-13	SB-2-(2-3)	17	1.11	9.66	9.32	96.0	
L3729-14	COMP-1-(0-1)	18	1.13	9.9	9.22	92.2	
L3729-15	COMP-2-(1-2)	19	1.18	9.69	9.12	93.3	
L3729-16	COMP-3-(2-3)	20	1.18	9.55	9.13	95.0	
L3729-17	COMP-4-(3-4)	21	1.14	9.62	8.92	91.7	
L3729-18	COMP-5-(4-5)	22	1.16	9.9	8.92	88.8	
L3731-01	POLY-DRUM	23	1.14	9.88	8.76	87.2	
L3740-01	SP-3A	24	1.18	9.92	8.57	84.6	
L3740-02	SP-3B	25	1.13	9.53	8.53	88.1	
L3740-03	SP-3C	26	1.17	9.59	7.63	76.7	
L3740-04	SP-3D	27	1.18	9.96	8.56	84.1	
L3742-01	TP-10-S1	28	1.18	9.6	8.73	89.7	
L3742-02	TP-10-S2	29	1.19	9.86	8.86	88.5	
L3742-03	TP-10-S	30	1.14	9.73	8.73	88.4	
L3742-05	TP-10-D1	31	1.1	9.69	8.41	85.1	
L3742-06	TP-10-D2	32	1.15	9.93	9.00	89.4	
L3742-07	TP-10-D	33	1.11	9.78	8.85	89.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$