



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/27/2020

OVENTEMP IN Celsius(°C): 106
Time IN: 13:00
In Date: 08/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110696

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3807-01	BFSN6	1	1.17	9.57	8.7	89.6	
L3807-02	BFSN7	2	1.17	9.71	7.75	77.0	
L3807-03	BFSN8	3	1.16	9.93	7.94	77.3	
L3807-04	BFSN9	4	1.17	9.51	8.04	82.4	
L3807-05	BFSP0	5	1.17	9.54	7.39	74.3	
L3807-06	BFWR1	6	1.16	9.63	8.43	85.8	
L3807-08	BFWS3	7	1.18	9.52	7.75	78.8	
L3807-09	VHBLK01	8	1.00	2.00	2.00	100.0	VHBLK
L3807-11	BFSSR8	9	1.17	9.6	8.82	90.7	
L3807-12	BFSSR9	10	1.15	9.52	8.79	91.3	
L3807-13	BFSS0	11	1.18	9.87	8.83	88.0	
L3807-14	BFSS1	12	1.19	9.91	8.57	84.6	
L3807-15	BFSS2	13	1.19	9.67	8.03	80.7	
L3807-16	BFWR2	14	1.17	9.74	9.00	91.4	
L3807-18	BFWS6	15	1.19	9.87	8.45	83.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$