



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/27/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:52
In Date: 08/26/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/27/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110705

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3498-01	HR-06-082620	1	1.16	9.65	9.1	93.5	
L3498-02	HR-06-082620	2	1.18	9.67	8.85	90.3	
L3498-03	HR-06-082620-1	3	1.13	9.54	8.84	91.7	
L3498-04	HR-06-082620-2	4	1.19	9.57	8.79	90.7	
L3803-01	KMH392P-END-1	5	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3803-02	KMH392P-END-2	6	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3803-03	BC260025COK-END-1	7	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3803-04	BC260025COK-END-2	8	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3804-01	GHA542P-END-1	9	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3804-02	GHA542P-END-2	10	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3804-03	HEH756Q-END-1	11	1.00	2.00	2.00	100.0	WOOD SAMPLE
L3804-04	HEH756Q-END-2	12	1.00	2.00	2.00	100.0	WOOD SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$