



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/31/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 08/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110734

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3508-25	TR-08-082820	44	1.15	9.94	9.6	96.1	
L3508-26	TR-08-082820	45	1.11	9.64	9.09	93.6	
L3508-27	TR-08-082820-1	46	1.19	9.95	9.57	95.7	
L3814-01	64-90-91	1	1.13	9.75	9.41	96.1	
L3814-05	31190	2	1.18	9.53	7.76	78.8	
L3814-06	31191	3	1.16	9.57	9.34	97.3	
L3814-07	31164	4	1.11	9.96	9.7	97.1	
L3830-01	CONCRETE-RO	47	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L3831-14	OF-3-20200827	5	1.19	9.85	5.97	55.2	
L3831-15	L3831-14MS	6	1.19	9.85	5.97	55.2	
L3831-16	L3831-14MSD	7	1.19	9.85	5.97	55.2	
L3831-17	OF-4-20200827	8	1.11	9.98	7.47	71.7	
L3831-18	DUP-02-20200827	9	1.13	9.84	6.3	59.4	
L3833-01	12018	41	1.18	9.71	9.22	94.3	
L3833-02	12018	42	1.18	9.59	9.13	94.5	
L3834-01	24745	43	1.14	9.79	5.43	49.6	
L3835-05	SVOC-GPC-BLANK	10	1.00	2.00	2.00	100.0	
L3835-06	PEST-GPC-BLANK	11	1.00	2.00	2.00	100.0	
L3835-07	PEST-GPC-BLANK-SPIKE	12	1.00	2.00	2.00	100.0	
L3835-08	PCB-GPC-BLANK	13	1.00	2.00	2.00	100.0	
L3835-09	PCB-GPC-BLANK-SPIKE	14	1.00	2.00	2.00	100.0	
L3835-10	SVOC-GPC2-BLANK	15	1.00	2.00	2.00	100.0	
L3835-11	PEST-GPC2-BLANK	16	1.00	2.00	2.00	100.0	
L3835-12	PEST-GPC2-BLANK-SPIKE	17	1.00	2.00	2.00	100.0	
L3835-13	PCB-GPC2-BLANK	18	1.00	2.00	2.00	100.0	
L3835-14	PCB-GPC2-BLANK-SPIKE	19	1.00	2.00	2.00	100.0	
L3836-01	PE-1	20	1.15	9.69	8.88	90.5	
L3836-02	PE-2	21	1.15	9.69	8.06	80.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/31/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 08/28/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/29/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110734

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3836-03	PE-3	22	1.12	9.56	8.15	83.3	
L3836-04	PE-4	23	1.14	9.86	8.89	88.9	
L3836-05	PE-5	24	1.18	9.69	8.81	89.7	
L3836-06	D82720	25	1.19	9.89	8.98	89.5	
L3837-06	789UMR-TP05-0612-01	26	1.13	9.85	9.38	94.6	
L3837-07	789UMR-TP05-1218-01	27	1.11	9.83	9.22	93.0	
L3837-08	789UMR-TP05-1824-01	28	1.16	9.79	9.27	94.0	
L3837-09	789UMR-TP05-3036-01	29	1.18	9.7	9.33	95.7	
L3837-10	789UMR-TP06-0006-01	30	1.19	9.63	8.46	86.1	
L3837-11	789UMR-TP07-1218-01	31	1.17	9.78	9.14	92.6	
L3837-12	789UMR-TP08-1824-01	32	1.18	9.6	9.00	92.9	
L3837-13	789UMR-TP11-0006-01	33	1.18	9.68	9.02	92.2	
L3837-14	L3837-13MS	34	1.18	9.68	9.02	92.2	
L3837-15	L3837-13MSD	35	1.18	9.68	9.02	92.2	
L3837-16	789UMR-TP11-0006-02	36	1.17	9.92	9.25	92.3	
L3837-17	789UMR-TP11-0612-01	37	1.15	9.72	9.3	95.1	
L3837-18	789UMR-TP11-2430-01	38	1.17	9.83	9.54	96.7	
L3837-19	789UMR-TP12-0002-01	39	1.14	9.8	9.47	96.2	
L3837-20	789UMR-TP13-0612-01	40	1.1	9.75	9.08	92.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$