



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 08/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 09/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110749

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3844-01	MC0AA3	1	1.11	9.91	4.63	40.0	
L3844-02	MC0AA4	2	1.16	9.94	4.18	34.4	
L3844-03	MC0AA5	3	1.12	9.58	4.74	42.8	
L3844-04	MC0AA6	4	1.13	9.9	4.48	38.2	
L3844-05	MC0AA7	5	1.15	9.56	5.17	47.8	
L3844-06	MC0AD2	6	1.17	9.71	4.75	41.9	
L3844-07	MC0AE0	7	1.18	9.83	4.29	36.0	
L3844-08	MC0AE3	8	1.15	9.95	4.81	41.6	
L3844-09	MC0AE5	9	1.16	9.89	4.75	41.1	
L3844-12	MC0AA8	10	1.19	9.92	6.42	59.9	
L3844-13	MC0AA9	11	1.16	9.8	5.86	54.4	
L3844-14	MC0AB0	12	1.13	9.82	7.19	69.7	
L3844-15	MC0AB1	13	1.13	9.53	6.89	68.6	
L3844-16	MC0AB2	14	1.16	9.98	6.69	62.7	
L3844-17	MC0AB3	15	1.15	9.94	4.93	43.0	
L3844-18	MC0AB5	16	1.18	9.96	4.23	34.7	
L3844-19	MC0AB7	17	1.12	9.77	4.99	44.7	
L3844-20	MC0AB8	18	1.13	9.83	5.74	53.0	
L3844-21	MC0AB8D	19	1.13	9.83	5.74	53.0	
L3844-22	MC0AB8S	20	1.13	9.83	5.74	53.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$