



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110753

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3506-11	OR-02-083120	80	1.18	9.66	9.02	92.5	
L3506-12	OR-02-083120	81	1.12	9.66	9.22	94.8	
L3506-13	OR-02-083120-1	82	1.1	9.66	8.59	87.5	
L3506-14	OR-02-083120-2	83	1.18	9.8	9.09	91.8	
L3828-02	BS-1	1	1.13	9.58	7.46	74.9	
L3828-03	BS-2	2	1.15	9.69	6.75	65.6	
L3828-04	BS-3	3	1.17	9.68	8.29	83.7	
L3846-05	BP-VPB181-GW-723-725	4	1.16	9.97	2.75	18.0	SLUDGE SAMPLE
L3847-01	LINDEN-JOP-BH-08-A	5	1.15	9.89	9.2	92.1	
L3847-02	LINDEN-JOP-BH-08-A	6	1.13	9.89	8.95	89.3	
L3847-03	LINDEN-JOP-BH-08-B	7	1.13	9.66	8.79	89.8	
L3847-04	LINDEN-JOP-BH-08-B	8	1.19	9.8	8.96	90.2	
L3847-05	LINDEN-JOP-BH-08-C	9	1.14	9.81	9.05	91.2	
L3847-06	LINDEN-JOP-BH-08-C	10	1.15	9.94	9.15	91.0	
L3847-07	LINDEN-JOP-BH-09-A	11	1.12	9.64	8.97	92.1	
L3847-08	LINDEN-JOP-BH-09-A	12	1.15	9.95	9.46	94.4	
L3847-09	LINDEN-JOP-BH-09-B	13	1.11	9.6	8.44	86.3	
L3847-10	LINDEN-JOP-BH-09-B	14	1.13	9.75	8.71	87.9	
L3847-11	LINDEN-JOP-BH-09-C	15	1.18	9.95	9.12	90.5	
L3847-12	LINDEN-JOP-BH-09-C	16	1.15	9.96	9.15	90.8	
L3847-13	LINDEN-JOP-BH-10-A	17	1.19	9.64	9.06	93.1	
L3847-14	LINDEN-JOP-BH-10-A	18	1.18	9.55	9.28	96.8	
L3847-15	LINDEN-JOP-BH-10-B	19	1.16	9.94	9.33	93.1	
L3847-16	LINDEN-JOP-BH-10-B	20	1.16	9.61	8.99	92.7	
L3847-17	LINDEN-JOP-BH-10-C	21	1.11	9.65	9.11	93.7	
L3847-18	LINDEN-JOP-BH-10-C	22	1.15	9.61	8.99	92.7	
L3847-19	LINDEN-JOP-BH-11-A	23	1.16	9.86	8.94	89.4	
L3847-20	LINDEN-JOP-BH-11-A	24	1.19	9.89	9.16	91.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110753

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3847-21	LINDEN-JOP-BH-11-B	25	1.16	9.64	8.53	86.9	
L3847-22	LINDEN-JOP-BH-11-B	26	1.13	9.53	8.47	87.4	
L3847-23	LINDEN-JOP-BH-11-C	27	1.15	9.97	9.19	91.2	
L3847-24	LINDEN-JOP-BH-11-C	28	1.18	9.98	9.36	93.0	
L3847-25	LINDEN-JOP-BH-12-A	29	1.13	9.61	8.75	89.9	
L3847-26	LINDEN-JOP-BH-12-A	30	1.15	9.59	8.76	90.2	
L3847-27	LINDEN-JOP-BH-12-B	31	1.17	9.63	8.67	88.7	
L3847-28	LINDEN-JOP-BH-12-B	32	1.16	9.59	8.62	88.5	
L3847-29	LINDEN-JOP-BH-12-C	33	1.14	9.66	8.7	88.7	
L3847-30	LINDEN-JOP-BH-12-C	34	1.13	9.71	8.69	88.1	
L3848-01	LINDEN-BH-03-C	35	1.11	9.74	9.06	92.1	
L3848-02	LINDEN-BH-03-C	36	1.16	9.69	9.12	93.3	
L3848-03	LINDEN-BH-02-A	37	1.16	9.75	8.81	89.1	
L3848-04	LINDEN-BH-02-A	38	1.15	9.96	9.22	91.6	
L3848-05	LINDEN-BH-02-B	39	1.15	9.96	9.55	95.3	
L3848-06	LINDEN-BH-02-B	40	1.16	9.77	9.19	93.3	
L3848-07	LINDEN-BH-02-C	41	1.15	9.94	9.09	90.3	
L3848-08	LINDEN-BH-02-C	42	1.15	9.69	8.82	89.8	
L3848-09	LINDEN-BH-01-A	43	1.19	9.72	9.06	92.3	
L3848-10	LINDEN-BH-01-A	44	1.1	9.8	9.16	92.6	
L3848-11	LINDEN-BH-04-A	45	1.18	9.77	9.25	93.9	
L3848-12	LINDEN-BH-04-A	46	1.18	9.97	9.51	94.8	
L3848-13	LINDEN-BH-04-B	47	1.14	9.91	9.31	93.2	
L3848-14	LINDEN-BH-04-B	48	1.15	9.66	9.04	92.7	
L3848-15	LINDEN-BH-04-C	49	1.1	9.57	8.53	87.7	
L3848-16	LINDEN-BH-04-C	50	1.19	9.67	8.66	88.1	
L3848-17	LINDEN-BH-03-A	51	1.17	9.82	9.29	93.9	
L3848-18	LINDEN-BH-03-A	52	1.14	9.88	9.39	94.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/1/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/31/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/01/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110753

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3848-19	LINDEN-BH-03-B	53	1.13	9.71	9.3	95.2	
L3848-20	LINDEN-BH-03-B	54	1.18	9.7	9.3	95.3	
L3848-21	LINDEN-BH-01-B	55	1.17	9.8	9.09	91.8	
L3848-22	LINDEN-BH-01-B	56	1.11	9.61	8.93	92.0	
L3848-23	LINDEN-BH-01-C	57	1.15	9.96	9.16	90.9	
L3848-24	LINDEN-BH-01-C	58	1.15	9.98	9.29	92.2	
L3848-25	LINDEN-FENCE-01	59	1.00	2.00	2.00	100.0	FENCE SAMPLE
L3848-26	LINDEN-FENCE-02	60	1.00	2.00	2.00	100.0	FENCE SAMPLE
L3854-01	DRUM-1-3	75	1.13	9.71	9.66	99.4	
L3854-02	DRUM-1	76	1.15	9.53	9.5	99.6	
L3854-03	DRUM-2	77	1.13	9.59	9.56	99.6	
L3854-04	DRUM-3	78	1.00	2.00	2.00	100.0	OILY DEBRIS
L3854-05	VNJ216-355	79	1.11	9.79	8.35	83.4	
L3857-01	PE-1A	61	1.16	9.83	8.9	89.3	
L3857-02	PE-2A	62	1.11	9.5	8.53	88.4	
L3857-03	PE-3A	63	1.16	9.62	8.7	89.1	
L3857-04	PE-4A	64	1.16	9.97	9.12	90.4	
L3857-05	PE-4B	65	1.14	9.94	9.37	93.5	
L3857-06	PE-1B	66	1.12	9.93	9.33	93.2	
L3857-07	PE-2B	67	1.1	9.59	9.12	94.5	
L3857-08	PE-3B	68	1.13	9.77	9.42	95.9	
L3857-09	PE-5	69	1.14	9.98	9.66	96.4	
L3857-10	OIL	70	1.00	2.00	2.00	100.0	OIL SAMPLE
L3857-14	DUP	71	1.13	9.86	9.00	90.1	
L3859-01	TP-1-S1	72	1.17	9.9	9.1	90.8	
L3859-02	TP-1-S2	73	1.16	9.69	8.88	90.5	
L3859-03	TP-1-S	74	1.17	9.57	8.8	90.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$