



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 09/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110800

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3746-08	359	1	1.19	9.72	6.2	58.7	
L3746-09	224	2	1.11	9.84	7.39	71.9	
L3871-01	TR-05-090220	3	1.13	9.64	9.08	93.4	
L3871-02	TR-05-090220	4	1.16	9.88	9.41	94.6	
L3871-03	TR-05-090220-1	5	1.1	9.81	9.31	94.3	
L3877-01	TP-1	6	1.13	9.68	8.6	87.4	
L3877-02	TP-1	7	1.16	9.76	8.01	79.7	
L3877-04	TP-1	8	1.16	9.76	8.01	79.7	
L3877-06	TP-2	9	1.16	9.98	8.21	79.9	
L3877-07	TP-2	10	1.13	9.81	8.93	89.9	
L3877-09	TP-2	11	1.13	9.81	8.93	89.9	
L3877-11	TP-3	12	1.15	9.85	8.97	89.9	
L3877-12	TP-3	13	1.13	9.59	8.22	83.8	
L3877-14	TP-3	14	1.13	9.59	8.22	83.8	
L3877-16	TP-4	15	1.14	9.97	9.06	89.7	
L3877-17	TP-4	16	1.16	9.94	8.82	87.2	
L3877-19	TP-4	17	1.16	9.94	8.82	87.2	
L3886-01	TOTE-31230	18	1.00	2.00	2.00	100.0	OIL SAMPLE
L3886-05	31227	19	1.00	2.00	2.00	100.0	OILY DEBRIS
L3886-06	31233	20	1.00	2.00	2.00	100.0	OILY-DEBRIS
L3886-07	31229	21	1.00	2.00	2.00	100.0	OILY DEBRIS
L3886-08	31228	22	1.00	2.00	2.00	100.0	OILY-DEBRIS
L3886-09	DRUM-1-COMP	23	1.00	2.00	2.00	100.0	OILY DEBRIS
L3887-01	CONCRETE-PILE	24	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L3887-02	CONCRETE-PILE	25	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L3888-01	9-2-TP-1	26	1.16	9.85	9.26	93.2	
L3888-02	9-2-TP-1	27	1.15	9.94	9.03	89.6	
L3888-04	9-2-TP-1-EPH-2	28	1.16	9.78	8.98	90.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/3/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 09/02/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/03/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110800

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3888-05	9-2-TP-2	29	1.13	9.79	8.57	85.9	
L3888-06	9-2-TP-2	30	1.14	9.91	8.78	87.1	
L3888-08	9-2-TP-2-EPH-2	31	1.13	9.61	8.44	86.2	
L3888-09	9-2-TP-3	32	1.1	9.6	8.78	90.4	
L3888-10	9-2-TP-3	33	1.13	9.94	9.06	90.0	
L3888-12	9-2-TP-3-EPH-2	34	1.18	9.83	8.84	88.6	
L3888-13	9-2-TP-4	35	1.11	9.98	9.3	92.3	
L3888-14	9-2-TP-4	36	1.14	9.92	8.99	89.4	
L3888-16	9-2-TP-4-EPH-2	37	1.13	9.88	9.11	91.2	
L3895-01	SB-6A- (7.5-8)	38	1.15	9.96	7.91	76.7	
L3895-02	SB-8A- (7-7.5)	39	1.13	9.78	7.42	72.7	
L3895-03	SB-9A- (7-7.5)	40	1.16	9.69	8.84	90.0	
L3895-04	SB-17W- (4.5-5)	41	1.18	9.68	8.75	89.1	
L3895-05	SB-17W- (7.5-8)	42	1.15	9.96	8.52	83.7	
L3895-06	SB-21A- (8.2-8.7)	43	1.19	9.74	7.67	75.8	
L3895-07	SB-22A- (8.5-9)	44	1.16	9.93	7.03	66.9	
L3895-08	SB-23A- (1.5-2)	45	1.16	9.76	6.79	65.5	
L3895-09	SB-23A- (7-7.5)	46	1.11	9.98	7.11	67.6	
L3895-10	SB-25A- (8.5-9)	47	1.15	9.96	8.12	79.1	
L3896-01	WC-1	48	1.13	9.59	7.02	69.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$