



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/14/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 09/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110835

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
L3808-04	PAV-B44-090320-0-1	1	1.18	9.68	9.38	96.5	
L3808-05	PAV-B44-090320-0-10	2	1.19	9.84	9.27	93.4	
L3808-06	PAV-B44-090320-19-20	3	1.14	9.82	9.13	92.1	
L3808-07	PAV-B44-090320-10-20	4	1.12	9.68	8.96	91.6	
L3808-08	PAV-B44-090320-DP	5	1.14	9.96	9.53	95.1	
L3808-09	PAV-B44-090320-DP	6	1.15	9.96	9.16	90.9	
L3826-07	BKG-1- (0.0-0.5)	39	1.13	9.55	7.67	77.7	
L3826-08	BKG-2- (0.0-0.5)	40	1.1	9.6	7.96	80.7	
L3826-09	BKG-3- (0.0-0.5)	41	1.18	9.73	7.92	78.8	
L3872-01	BU-03-090420	42	1.13	9.8	8.14	80.9	
L3872-02	BU-03-090420	43	1.18	9.51	7.00	69.9	
L3872-03	BU-03-090420-1	44	1.16	9.93	9.18	91.4	
L3912-02	LS-1-B-2-ELEVATOR-SHAFT	7	1.00	2.00	2.00	100.0	OIL SAMPLE
L3915-01	WC-PM-A1-0-4-C	8	1.18	9.98	8.14	79.1	
L3915-05	WC-PM-E-0-5-C	9	1.14	9.9	8.7	86.3	
L3915-09	WC-PM-B2-0-3.5-C	10	1.13	9.87	9.24	92.8	
L3915-13	WC-PM-A2-10-14-C	11	1.14	9.65	8.23	83.3	
L3915-18	WC-PM-A2-10-14-D	12	1.12	9.57	8.19	83.7	
L3916-01	WC-1	13	1.16	9.55	8.33	85.5	
L3916-04	WC-2	14	1.18	9.96	8.96	88.6	
L3918-01	72-11986	45	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
L3919-01	TP-10	15	1.13	9.59	8.69	89.4	
L3919-02	TP-10	16	1.16	9.59	8.64	88.7	
L3919-04	TP-11	17	1.12	9.55	8.68	89.7	
L3919-05	TP-11	18	1.1	9.8	8.8	88.5	
L3919-07	TP-12	19	1.13	9.66	9.1	93.4	
L3919-08	TP-12	20	1.18	9.62	8.76	89.8	
L3919-10	TP-13	21	1.15	9.96	8.75	86.3	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/14/2020

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/04/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 09/05/2020
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB110835

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
L3919-11	TP-13	22	1.13	9.64	8.61	87.9	
L3919-13	TP-14	23	1.14	9.61	8.41	85.8	
L3919-14	TP-14	24	1.17	9.57	8.36	85.6	
L3919-16	TP-15	25	1.13	9.74	8.74	88.4	
L3919-17	TP-15	26	1.13	9.85	9.07	91.1	
L3923-01	BP-VPB181-GW-858-860	27	1.14	9.97	3.03	21.4	SLUDGE SAMPLE
L3923-07	BP-VPB181-GW-878-880	28	1.19	9.59	2.18	11.8	SLUDGE SAMPLE
L3925-05	SVOC-GPC-BLANK	29	1.00	2.00	2.00	100.0	
L3925-06	PEST-GPC-BLANK	30	1.00	2.00	2.00	100.0	
L3925-07	PEST-GPC-BLANK-SPIKE	31	1.00	2.00	2.00	100.0	
L3925-08	PCB-GPC-BLANK	32	1.00	2.00	2.00	100.0	
L3925-09	PCB-GPC-BLANK-SPIKE	33	1.00	2.00	2.00	100.0	
L3925-10	SVOC-GPC2-BLANK	34	1.00	2.00	2.00	100.0	
L3925-11	PEST-GPC2-BLANK	35	1.00	2.00	2.00	100.0	
L3925-12	PEST-GPC2-BLANK-SPIKE	36	1.00	2.00	2.00	100.0	
L3925-13	PCB-GCP2-BLANK	37	1.00	2.00	2.00	100.0	
L3925-14	PCB-GCP2-BLANK-SPIKE	38	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$